

PHILIPS

Disposable 30-day
fine filter

DreamStation Go

Filters dust, pollen, & more

Replace every 30-days

95% efficiency at $< .7$ micron

6-pack

HH1454



30-day fine filter

Philips 30-day disposable replacement fine filters help clean the air you breathe by capturing dust, pollen, mold spores, pet dander and more with a 95% efficiency at .5-.7 micron.

Specifications

Technical information

Replacement: Every 30 days
Filter material: Blended Synthetic Fiber
Filter efficiency: 95% Efficient @ .5-.7 micron

Packaging details

Content: Includes 2 filters

Part Number

Disposable 30-day fine filter: 1133744

Service

Warranty: 15-days

Compatibility

PAP System: DreamStation Go CPAP, AutoPAP

