

PHILIPS



Under the nose nasal mask

DreamWear Nasal

Under-the-nose nasal mask

4 different cushion sizes

Includes frame and headgear



HH1111/01

DreamWear

Nasal mask with mask cushion underneath the nose

The innovative design of DreamWear offers you the freedom and comfort to sleep in the way you would like.

Freedom of movement

- For a more natural sleep
- Soft, flexible mask frame
- Innovative mask design for more freedom in sleep therapy
- More comfort with DreamWear

Innovative design

- Stability like never before

No red marks on the bridge of the nose

- Cushion under the nose

Unobstructed view

- Provides freedom of movement throughout the night
- Open field of vision

Highlights

Nasal cushion




The mask cushion lies gently under the nose, which helps to prevent pressure points on the bridge of the nose and unpleasant skin irritations on the nostrils.

A more natural sleep



DreamWear combines the advantages of minimal contact masks and nasal masks, provides high comfort and can contribute to successful sleep therapy.

Top-of-head tube connection




The exhaled air is not discharged directly at the nose, but at the top of the head. This hose positioning allows for more freedom of movement throughout the night.

Read or watch TV



Thanks to the open design, DreamWear offers an open field of vision, so you can easily wear glasses and read, watch TV or use a tablet PC in bed before going to sleep.

Soft, flexible mask frame



The flexible material, through which the air flows, gently rests on the skin, prevents pressure points and provides a comfortable fit.

Innovative mask design

The innovative design of DreamWear helps to prevent pressure sores and skin irritation on the nostrils and on the nose bridge, and allows the user, in addition to a clear field of vision, more freedom of movement through the hose connection on the head.

More comfort with DreamWear

Users said DreamWear was more comfortable to wear than their previously prescribed mask, as it is barely noticeable.

Specifications

Cleaning

Hand wash with mild detergent: Rinse thoroughly, hang or air dry

Product Details

Materials: Silicone, Polycarbonate, Polyurethane foam, Spandex, Thermoplastic elastomer

Materials not included: Not made with natural rubber latex

Length of Medium Frame: 20 inches

Package contents: Medium frame, S, M, L, and MW Cushions, Headgear w/ arms, Fabric Wraps

Innovative design

Tube connection on top of head: Gives freedom of movement

Freedom of movement

Clear field of vision: Individual sleep position

Design

Easy to use: Easy assembly and disassembly

Comfort: Reduces nasal irritation

Stability: Back panel with crown strap

Stays in place: Non-slip headgear

Quick release: 360° rotating elbow

Cleaning: Fabric wraps

Hand wash with mild detergent: Rinse thoroughly, hang to dry¹

Service

Warranty: 90-days

Replacement

Inspect daily for wear: Replace at first sign of wear



* *Source: Philips User Survey 2015

* The costs of DreamWear are private and are usually not reimbursed by the statutory health insurance.

* The scope of supply includes a standard mask frame (size M) and four cushions for all sizes (S, M, MW, L). The sleep therapy device is not included.