

PHILIPS

Efficia

Patient monitors



Efficia patient monitor product structure



Model	Part number	Screen	CM10 863301 10"	CM12 863303 12"	CM100 863300 10"	CM120 863302 12"	CM150 863304 15"
Measurement parameter options (accessories must be ordered separately)							
3- and 5-lead ECG			•	•	•	•	•
12-lead ECG	B07			□		□	□
Philips ST analysis			•	•	•	•	•
Impedance respiration			•	•	•	•	•
Philips SpO ₂ ¹	B01		□	□	□	□	□
Masimo SET ²	B02		□	□	□	□	□
Nellcor OxiMax SpO ₂ ²	B11		□	□	□	□	□
Masimo rainbow SpHb/SpOC ³	B03				□	□	□
Masimo rainbow SpCO ³	B04				□	□	□
Masimo rainbow RRa ³	B05				□	□	□
Philips NBP			•	•	•	•	•
Dual continuous temp			•	•	•	•	•
Dual IBP including overlapping IBP display software	B10		□	□	□	□	□
CO ₂ Oridion Microstream ⁴	B20		□	□	□	□	□
CO ₂ Respironics Capnostat Mainstream ⁴	B21		□	□	□	□	□
CO ₂ Respironics CapnoTrak Sidestream ⁴	B22		□	□	□	□	□
Cardiac output ⁵ (with hemodynamics calculations and wedge pressure)	B30			□		□	□

• Standard □ Optional

Model		CM10	CM12	CM100	CM120	CM150
Software options						
Philips basic arrhythmia		•	•	•	•	•
Philips enhanced arrhythmia	C01	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Assisting venous puncture	C04	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Barcode reader support ⁶	C08	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Philips ST Map	C10		<input type="checkbox"/>		•	•
48 hours of full disclosure (all waveforms)	C15	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
HL7 output	C20	<input type="checkbox"/>	<input type="checkbox"/>	•	•	•
HL7 inbound ⁷	C21			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Features and hardware options						
10-inch display		•		•		
12-inch display			•		•	
15-inch display						•
Touchscreen	E01	<input type="checkbox"/>	<input type="checkbox"/>	•	•	•
Maximum number of waveforms displayed		8	12	8	8	12
External display connection		•	•	•	•	•
Internal recorder ⁸	E10	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Display mini trend, 8 hours			•		•	•
USB ports		2	2	2	2	2
Night mode		•	•	•	•	•
ECG analog output	E03	<input type="checkbox"/>	•	•	•	•
Dual speaker				•	•	•
HL7 wireless LAN connectivity	E20	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Dashboard connectivity		•	•	•	•	•
Bed rail hook	E16	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Accessories and cable options						
Accessories not included	A02	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Adult accessories kit (3-lead) ^{9,10}	A01	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pediatric accessories kit (3-lead) ^{10,11}	A03	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Neonatal accessories kits (3-lead) ^{10,12}	A04	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Adult accessories kit (5-lead) ^{9,10}	A07	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Adult premium accessories kit (3-lead) ^{2,10}	A05	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Adult premium accessories kit (5-lead) ^{2,10}	A06	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CO ₂ Mainstream sensor ¹³	A21	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Reusable adult/pediatric CO ₂ Mainstream adapter ¹⁴	K21	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Battery options						
3-cell lithium-ion battery	E11	<input type="checkbox"/>	<input type="checkbox"/>			
9-cell lithium-ion battery	E12	<input type="checkbox"/>	<input type="checkbox"/>	•	<input type="checkbox"/>	<input type="checkbox"/>
Dual 3-cell lithium-ion batteries	E13		<input type="checkbox"/> ¹⁵			
Dual 9-cell lithium-ion batteries	E14		<input type="checkbox"/> ¹⁵		<input type="checkbox"/>	<input type="checkbox"/>

• Standard Optional

1. Default unless Masimo or Nellcor selected.
2. May not be available in all countries.
3. Can only be ordered in combination with B02.
4. Only one of B20, B21 or B22 can be ordered.
5. Can only be ordered in combination with B10.
6. Barcode reader must be ordered separately:
- 98980 317 6611
- 98980 314 7821
7. Philips IntelliBridge Enterprise Interoperability solution required.

8. Can only be ordered together with either option E12 or option E14.
9. 3-lead ECG one-piece trunk cable (grabber), adult SpO₂ sensor, adult NBP cuff and NBP hose.
10. SpO₂ sensor and cable only included when B01 is selected.
11. 3-lead ECG one-piece trunk cable (grabber), pediatric SpO₂ sensor, pediatric NBP cuff and NBP hose.

12. 3-lead trunk cable, disposable neonatal ECG lead wires/electrodes, Philips neonatal hand/foot SpO₂ sensor, SpO₂ extension cable, neonatal NBP hose, 25 (5 cuffs per each size) neonatal single-patient NBP cuffs (various sizes).
13. Can only be ordered in combination with B21 Respiration Capnostat Mainstream sensor M2501A.
14. Can only be ordered in combination with A21.
15. Can only be ordered in combination with E01.

The following are trademarks or registered trademarks of Masimo Corporation: DCI, FastSat, LNCS, LNOP, Masimo, PVI, SET, SpCO, SpHb, SpMet, TF-I, rainbow ReSposable, APOD, rainbow, M-LNCS, rainbow Acoustic Monitoring, rainbow Acoustic Sensor, rainbow Acoustic.

Oridion is a trademark and Microstream is a registered trademark of Oridion Medical Ltd.

Nellcor and OxiMax are trademarks of Medtronic.

© 2018 Koninklijke Philips N.V. All rights are reserved. Philips reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.



philips.com

Printed in The Netherlands.
4522 991 35541 * JUN 2018