

# PHILIPS



Headgear (Standard)

Amara Gel

Standard size



HH1220

## When it's right, it just clicks

The Amara full face mask has significantly fewer parts than leading traditional full-face masks. It requires just a single click to disassemble and reassemble the headgear.

### Minimalistic design

- 60% fewer parts

## Highlights

### Simplicity

With only four parts, assembly and cleaning are quick and easy.

## Specifications

### Specifications

**Headgear material:** UBL/Urethane foam/nylon and Lycra

**Cleaning: headgear, tubing:** Hand wash with mild detergent, Rinse thoroughly, hang to dry

### Design

**Comfort:** Self-adjusting

**Fewer parts:** Comfort

### Consumer Trade Item

**GTIN:** 606959007185

**Country of origin:** CN

**Harmonized system code:** 90189084

### Outer Carton

**GTIN:** 20075020086925

### Service

**Warranty:** 90-days

### Other

**GTIN:** 50075020098202

