

# PHILIPS

IPL hair removal  
system

Lumea Comfort

For use on body and face

22 minutes to treat lower legs

Lifetime >100.000 light  
pulses

SC1981/00



## Compact IPL solution

to prevent hair regrowth for smooth skin

Philips Lumea Comfort IPL hair removal system works wonders to prevent the reappearance of hair on face as well as body. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

### Long lasting smoothness

- Effortlessly effective
- Enjoy smooth skin every day

### Hair removal with IPL technology

- Adapted for safe use in the comfort of your own home

### Most compact Lumea design

- No replacement parts, no hidden costs
- The smallest and lightest device within the Lumea range

### Safe and effective on face and body

- Five adjustable light energy settings
- Safe and effective even on sensitive areas
- Skin Color Sensor
- Slide and Flash mode for a more convenient application
- Suitable for different hair and skin types


# Highlights

## Enjoy smooth skin every day



Philips Lumea applies gentle pulses of light to the hair root. This results in the hair shedding naturally and hair regrowth is inhibited. Repeating this treatment regularly will keep your skin smooth everyday.

## Effortlessly effective



Our clinical studies have shown significant hair reduction in as little as four bi-weekly treatments resulting with smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual hair regrowth. - To enlarge this picture, please click on the image in the image gallery on the top of this page

## IPL technology



Philips Lumea Comfort uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for effective use in the safety and comfort of your own home. Philips worked closely with skin care experts to develop this technology. For more than 10

years Philips carried out extensive consumer research on IPL technology with more than 2000 women.

## Suitable hair and skin types

Haarfarbe Haarkleur Hair color Couleur de poils	Hautton / Huidskleur / Skin tone / Couleur de la peau					
	Light	Light-Medium	Medium	Medium-Dark	Dark	Very Dark
Black	✓	✓	✓	✓	✗	✗
Dark Brown	✓	✓	✓	✓	✗	✗
Medium Brown	✓	✓	✓	✓	✗	✗
Light Brown	✓	✓	✓	✓	✗	✗
Red	✗	✗	✗	✗	✗	✗
Light-Blonde	✗	✗	✗	✗	✗	✗
White/Grey	✗	✗	✗	✗	✗	✗

Philips Lumea Comfort works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Lumea Comfort is NOT effective on red, light-blonde or white/grey hair. - To enlarge this picture, please click on the image in the image gallery on the top of this page

## Application areas



Use with confidence to remove hair on face (upper lip, chin ) and body parts including legs, arm pits, bikini area, belly and arms.

## 5 light energy settings




Philips Lumea Comfort has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is safe and gentle to use even on sensitive skin and sensitive body areas.

## Skin Color Sensor



The unique Skin Color Sensor measure the treated skin complexion at the beginning of each session and occasionally during the session. If it detects a skin tone that is too dark for Lumea Comfort, it will automatically stop emitting pulses.

## Slide and Flash mode



Lumea Comfort enables a sliding motion for a more convenient and continuous application.

## Smallest and lightest Lumea

The Philips Lumea Comfort is the smallest and lightest in the Lumea range. It is easy to hold and easy to reach all the areas you want treat.

## No replacement parts needed



Philips Lumea Comfort is a complete solution that comes ready to work straight out of the box. More importantly, it does not require any replacement lamps or gels.

# Specifications

## Power

Voltage: 100-240 V

## Service

Warranty: 2 years global warranty

## Technical specifications

**High performance lamp:** Delivers <gt/> 100,000 flashes

## Application areas

**Body areas:** Arms, Legs, Belly, Bikini, Armpits

**Face areas:** Chin, Sideburns, Upper lip

## Items included

**Instructions for use:** Quick Start Guide, User manual

**Adapter:** 12V / 1500mA

## Safety and adjustable settings

**Integrated UV filter:** Protects skin from UV light

**5 light energy settings:** Adjustable to your skin type

**Integrated safety system:** Prevent unintentional flashing

**Skin Tone sensor:** Detects your skin tone

## Application time

**Lower leg:** 11 min

**Armpit:** 1,5 min

**Bikini line:** 1,5 min

**Face areas:** 2,5 min

## Application mode

**Slide and Flash:** For easy motion application

**Step and Flash:** For treatment on curvy areas

**Corded / cordless use:** Corded

