

# PHILIPS

Under-the-nose  
mask

DreamWear Nasal

Top-of-head tube design

Enhanced fit and seal

Updated headgear with arms  
design

Available in multiple sizes



HH1111/01



## Like wearing nothing at all\*

Nasal mask with mask cushion underneath the nose

Philips Respironics DreamWear, now with updated headgear, delivers a new level of stability & performance to the mask. The DreamWear headgear with arm 's unique slip-resistant design helps keep the mask in place throughout the night.

#### **Innovative design**

- Enhanced fit & seal
- Sleep the way you dream about
- Stability like never before

#### **Modular design**

- Exceptional comfort

# Highlights

## Exceptional comfort

94% of DreamWear nasal users surveyed said they slept comfortably while wearing the DreamWear mask with updated headgear.<sup>1</sup> Users find DreamWear Nasal mask with updated headgear more comfortable than the ResMed Airfit N30i nasal mask.<sup>1</sup>

## Enhanced fit & seal

DreamWear nasal users<sup>1</sup> said that they were more satisfied with the fit & seal\*\* of the DreamWear mask with the updated headgear than their prescribed mask. According to DreamWear nasal users the DreamWear with updated headgear has a better fit than ResMed Airfit N30i Nasal mask.<sup>2</sup>

## Sleep the way you dream about

Keeping your mask & headgear in place means a higher quality of sleep all night long. A majority of DreamWear nasal users<sup>1</sup>

surveyed stated that they were more satisfied with the quality of their sleep with the DreamWear mask with updated headgear. Unrestricted by cumbersome tubing, patients can sleep in any position<sup>1</sup> they desire – on their stomach, side, or back.

## Stability like never before

With an innovative design that helps to reduce slippage, DreamWear headgear with arms delivers additional stability. The arms work the same way as arms on eye glasses by keeping the mask & headgear from moving when being used. 88% of DreamWear nasal users<sup>1</sup> said our mask with the updated headgear was more stable than the original headgear.

# Specifications

## Cleaning

**Hand wash with mild detergent:** Rinse thoroughly, hang or air dry

## Product Details

**Materials:** Silicone, Polycarbonate, Polyurethane foam, Spandex, Thermoplastic elastomer

**Materials not included:** Not made with natural rubber latex

**Length of Small Frame:** 19 inches

**Length of Medium Frame:** 20 inches

**Package contents:** Medium frame, S, M, L, and MW Cushions, Headgear w/ arms, Fabric Wraps

## Service

**Warranty:** 90-days

## Warning: Magnetic Interference

**Beware with pacemakers:** Consult your physician before use.



\* 1. Patient Preference study with DreamWear Nasal with original headgear users (n=140). Jan/March 2019 at 30 days. Original headgear references to the slit headgear design on DreamWear UTN Nasal mask  
\* 2. A 30-day Patient Preference study with DreamWear UTN Nasal original headgear users; n = 121 conducted April / May 2019  
\* \*\*Fit and seal were surveyed independent of each other.  
\* \*2015 Philips User Preference Questionnaire