

PHILIPS

IPL - Hair removal
device

Lumea Prestige

For use on body and face

15 minutes to treat lower legs

>250,000 lamp flashes

SC2007/00




Highlights

IPL technology



Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 volunteers.

Enjoy smooth skin every day



Philips Lumea applies gentle pulses of light to the hair root. As a consequence, the amount of hair your body grows gradually decreases. Repeating the treatment leaves your skin beautifully hair free and touchably smooth.

Effortlessly effective

Clinical results after treatment: **example 1** (leg)



Clinical results after treatment: **example 2** (bikini line)



Our clinical studies have shown significant hair reduction in as little as four fortnightly treatments resulting in smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between

treatments may vary based on your individual hair regrowth. - To enlarge this picture, please click on the image in the image gallery on the top of this page

Suitable hair and skin types




Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blonde or white/grey hair. Philips Lumea is also NOT suitable for very dark skin.

Application areas



Use with confidence to remove hair on face (upper lip, chin and sideburns) and body parts including legs, underarm, bikini area, belly and arms.

Fast treatment on legs



Big treatment window for fast application on large areas like legs.

Extra safety for the face




The precision attachment is equipped with an extra integrated light filter and enables a safe and convenient application on the upper lip, chin and sideburns. It can also be used on other sensitive areas on the body, such as the underarms, the bikini area and hard-to-reach areas.

Gentle treatment



Gentle to your skin.

Cordless convenience



Lumea's unique cordless design is perfect for maximum reach and ease of use.

Specifications

Power

Voltage: 100-240 V

Service

Warranty: 2 years global warranty

Technical specifications

High performance lamp: Delivers <gt/> 250,000 flashes

Application areas

Body areas: Arms, Legs, Belly, Bikini, Armpits

Face areas: Chin, Sideburns, Upper lip

Attachments

Precision attachment (2cm²): For safe application on face

Body attachment (4 cm²): For use below the neckline

Items included

Instructions for use: User manual

Charger: Multi voltage charger

Storage: Luxurious pouch

Accessory: Cleaning cloth

Power

Charging time: 100 min

Running time: Minimum 370 flashes

Safety and adjustable settings

Integrated UV filter: Protects skin from UV light
5 light energy settings: Adjustable to your skin type

Integrated safety system: Prevent unintentional flashing

Application time

Lower leg: 8 min

Armpit: 1 min

Bikini line: 1 min

Face areas: 1 min

Application mode

Slide and Flash: For convenient application

Corded/cordless use: Cordless

Stamp and Flash: For treatment on curvy areas

