

PHILIPS

IMT


Threshold

Inspiratory Muscle Trainer

Increases endurance

Increases respiratory strength

HH1332



Pressure that's good for you

With a flow-dependent, one-way valve, threshold devices maintain a constant pressure, regardless of the patient's airflow.

Convenient design

- Reliable and easy to clean
- Works in any position

Mouthpiece/mask option

- Mouthpiece/mask option

Highlights

Mouthpiece/mask option

Mouthpiece/mask option for flexibility

Reliable and easy to clean

Threshold is made of durable, high-impact acrylic.

Works in any position

Effective therapy in any position

Specifications

Product details

Warranty: Replace after 1 year

Intended Use: As directed by your physician

Cleaning: Wash in warm soapy water after each use. Flush all parts with clean water.

Shake out excess water and air dry. **DO NOT BOIL OR HEAT.**

Includes: training diary

