

# PHILIPS

Compressor  
Nebulizer System

InnoSpire Deluxe

Robust Home use device

Treatment in under 5 minutes\*

\*Using a 2.5ml Salbutamol dose

HH1350/00



## Reliable therapy with built-in convenience

InnoSpire Deluxe is an easy to use and reliable compressor nebulizer system that includes built-in storage to keep all your therapy accessories and supplies organized and handy.

### For use in hospital and home

- Designed for continuous use

### Short treatment times

- Treatment delivered in under 5 minutes\*

### Designed for convenience

- Sturdy handle for easy transportation
- Everything you need in one place

### For commonly prescribed medications

- Compatible with commonly prescribed nebulizer medications

# Highlights

## Convenience meets reliability



InnoSpire Deluxe is designed to make your nebuliser treatment convenient and easy, with added built-in storage to house all your accessories.

## For home use



Robust enough for continuous use in both a clinical and family environment, InnoSpire Elegance provides consistent and reliable nebulizer therapy.

## Robust Design




Weighing just 4.4 lbs, InnoSpire Deluxe is light enough to pick up and move from room to room at home. With your supplies stored neatly away within the device, your nebulizer equipment stays safely together in one place.

## Treatment time



Combined with our reusable SideStream nebulizer, it is designed to provide fast drug delivery and short treatment times.

## Versatility



Capable of nebulizing a wide variety of medications that are routinely used to treat various respiratory conditions.

# Specifications

## Thermal overload protector Warranty

**Service**  
Warranty: 3 years

**Design**  
Convenient storage compartment: 2.0 kg (4.4 lbs)

**Operating parameters**  
Maximum flow rate: 10.78 LPM  
Maximum pressure: 268 kPa (35 psi)  
Electricity supply: 120V/60Hz, 1.7 Amp  
Flow rate: 7 LPM @ 70 kPa (0.7 bar)

**Product details**  
Weight: 2.2 kg  
Size: 380 x 170 x 110 mm  
Noise level: 58 ±3 dBA  
Treatment time (under 5 min\*): \* 2.5ml Salbutamol dose

