

PHILIPS

Reusable Nebulizer

SideStream

Treatment in 6 minutes*

Can be used with mask

Suitable for adults & children

HH1320/00



Highlights

Effective delivery over time



Staying compliant with your breathing treatment is easier when taking your medication is consistent and reliable. SideStream nebulizers keep the driving air and the liquid drug separate for ongoing performance.

Fast, effective treatment*



The nebulizer top, combined with an InnoSpire compressor, SideStream delivers faster aerosol drug delivery and shorter treatment times, so you can get back to your daily life.

Fast, effective treatment time



The SideStream nebulizer's venturi system boosts the airflow supplied by the compressor, delivering medication more efficiently and making your overall treatment faster.

Specifications

Maintenance

Replace: After 6 months

Product details

Aerosol flow to patient: 16 L/min

Driving gas: 6 L/min

Venturi air entrainment: 10 L/min

