

# PHILIPS

## IMT

### Threshold


Inspiratory Muscle Trainer  
device

Increases endurance

Increases respiratory strength

Asthma and COPD

HH1332



## Pressure that's good for you

With a flow-dependent, one-way valve, threshold devices maintain a constant pressure, regardless of the patient's airflow.

#### Convenient design

- Reliable and easy to clean
- Works in any position

#### Mouthpiece/mask option

- Mouthpiece/mask option

#### Intended use

- Intended use of IMT

## Highlights

**Mouthpiece/mask option**

Mouthpiece/mask option for flexibility

**Reliable and easy to clean**

Threshold is made of durable, high-impact acrylic.

**Works in any position**

Effective therapy in any position

**Intended use**

Threshold IMT inspiratory muscle trainer can be used for muscle strength and endurance.

## Specifications

**Product details**

**Warranty:** Replace after 1 year

**Intended Use:** As directed by your physician

**Cleaning:** Wash in warm soapy water after each

use. Flush all parts with clean water. Shake out excess water and air dry. DO NOT BOIL OR HEAT.

**Includes:** training diary

