



The Philips Mask Adhesive, a new accessory for AF531 and PerformaTrak oro-nasal mask securement

The Philips Mask Adhesive* is an accessory intended for single-patient use during the COVID-19 pandemic in the hospital environment to help secure oro-nasal mask cushions to a patient's skin during noninvasive ventilation.

As the clinical needs of COVID-19 patients evolve, Philips is steadfast in its commitment to provide support to health systems and care providers as they seek to optimize noninvasive ventilation treatment.

The Mask Adhesive is intended to adhere the mask to the patient's face, with the goal of reducing the force needed from the headgear to secure the mask.

Features and benefits



Overview

The Mask Adhesive seeks to help secure the mask cushion to the patient's face, and address clinicians' concern of mask overtightening, which can lead to skin injury over time.

It also allows for mask removal while the adhesive remains on the face for re-application, enabling the patient to periodically take medication or to eat.



Easy to use

It's easy to apply the Mask Adhesive using the mask forehead pad alignment feature, and it's just as easy to remove when you're ready.



Patient friendly

Created with patient comfort in mind, the Mask Adhesive can remain in place, while the mask is lifted, allowing the patient to eat, drink and take oral medication. The mask is then reapplied directly back onto the adhesive on the patient's face to continue noninvasive ventilation treatment.





Secure mask fit

Using both the Mask Adhesive and the headgear can help provide a more secure mask fit. The Mask Adhesive is designed to meet a specific need of providers treating COVID-19 patients to secure the mask to the patient's face.

The Mask Adhesive, with the AF531 or the PerformaTrak oro-nasal masks, can help clinicians deliver exceptional noninvasive ventilation. Using the mask adhesive can give clinicians confidence to explore all noninvasive ventilation therapy options to support high-quality patient care.

Details

- Medical grade acrylic adhesive, used for First Aid and wound care dressings
- Hypoallergenic and latex free
- Fluid resistant
- Conformable to the patient's facial structure
- Adhesive can be re-used until it no longer adheres, at which point it should be replaced
- Replace Mask Adhesive as needed (at least every 3 days)
- Mask Adhesive can be used with both CapStrap and 4 Point Headgear



Mask Adhesive on AF531



Mask Adhesive on PerformaTrak

Ordering Information

Item	Part number – 12NC	Part number – PRD	Mask compatibility
S, Mask Adhesive, 20pk	989805664451	1146923	AF531 / PerformaTrak
M, Mask Adhesive, 20pk	989805664461	1146924	AF531 / PerformaTrak
L, Mask Adhesive, 20pk	989805664471	1146925	AF531 / PerformaTrak

*The Philips Mask Adhesive product is being provided in the USA under the FDA Enforcement Policy for Ventilators During the COVID-19 Public Health Emergency Guidance. As such, it is for use only during the U.S. Health and Human Services COVID-19 public health emergency declaration, unless terminated or revoked earlier, after which the product may no longer be used.

© 2020 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.

www.philips.com/respironics



Caution: US federal law restricts these devices to sale by or on the order of a physician.

edoc RB 7/2/20 MCI 4109953 v00
1010 Murry Ridge Lane, Murrysville, PA 15668
800 345 6443 • 724 387 4000