




Philips AVENT
Manual breast pump

with VIA Storage System



SCF290/12

Express more. Store easily.

Express, store and feed in easy steps

A manual breast pump offers you total control of the suction level and pumping rhythm. The unique Philips AVENT Manual Breast Pump assures you enhanced efficiency, being clinically proven to express more milk than a hospital electric pump.

For safe storage and transport

- Leak-proof twist-on system

Fits all Philips AVENT breast pumps and teats

- Use adaptor to express and feed

Clinically proven effective

- Patented Let-down Massage Cushion

Stimulate and Express

- Soft massage cushions trigger natural let-down, like baby

Benefits and features:

- Express and store in one easy step
- Clinically proven

PHILIPS
AVENT

sense **and** simplicity

Specifications

What is included

- Manual breast pump: 1 pcs
- Pre-sterilized VIA Cup (180 ml/ 6 oz): 5 pcs
- Milk storage cup lid: 5 pcs
- Milk storage cup adaptor: 2 pcs
- Extra soft Newborn Flow nipple: 2 pcs
- Newborn nipple travel pack: 2 pcs
- Spare parts: 1 pcs

Country of origin

- England

Development stages

- Stages: 0 - 6 months

Highlights

Soft massage cushions



Patented massage cushions of the Philips AVENT breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

Leak-proof twist-on system




For safe storage and transport

Adaptor to express and feed



Use the adaptor to express and feed with the same cup.

Patented Let-down Massage Cushion



Five, soft petal massagers of the Philips AVENT breast pump gently flex in and out as you pump and work with the vacuum to imitate your baby's suckling.



Issue date 2013-02-15

Version: 7.3.1

12 NC: 0000 000 00000

© 2013 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com