

Philips AVENT
4-in-1 Electric Steam
Sterilizer

Electric Sterilizer

SCF287/02



The easiest way to sterilize

Adjustable size saves kitchen space

Our new electric steam sterilizer has been designed to make sterilizing bottles as simple as possible. The 4-in-1 sterilizer has an integrated dishwasher basket for hassle free pre-cleaning and an advanced display to keep you informed.

4-in-1 sterilizer

- Adjustable size and integrated dishwasher basket

Keeps you informed

- Advanced display with progress and cool-down indication

Fits all bottle types

- Suitable for wide neck and narrow neck bottles

Extra speed and safety

- 6 minute cycle with automatic shut-off

Effective sterilization

- Natural steam sterilization kills 99.9% of harmful germs

PHILIPS
AVENT

sense and simplicity

Highlights

4-in-1 sterilizer



Due to its adjustable size, this 4-in-1 Philips AVENT sterilizer perfectly fits the items you would like to sterilize and always takes up the least amount of kitchen space. The 4-in-1 functionality gives you: 1) Small size configuration for sterilizing soothers 2) Medium size configuration for sterilizing several bottles, toddler plates, knives, and forks 3) Large size configuration for sterilizing a full load of up to 6 bottles 4) Integrated dishwasher basket for hassle free pre-cleaning of baby bottles and accessories

Advanced display



The advanced display keeps you informed throughout the sterilization cycle and lets you know when the bottles are cooled down.


Fits all bottle types



The sterilizer has been designed to fit both wide neck and narrow neck bottles as well as different breast pumps and accessories. This gives you the freedom to choose the bottles and accessories that best suit you and your baby. At full capacity the

sterilizer holds six 11oz / 330ml Philips AVENT bottles.

6 minute cycle & auto shut-off



Your bottles and accessories will be sterilized in just 6 minutes. Once sterilization is complete the device will shut off automatically, ensuring extra safety and less energy consumption.

Kills 99.9% of harmful germs



Sterilizing protects your baby from particularly harmful bacteria until his immune system is strong enough. Philips AVENT Sterilizers use the hospital method of steam sterilization, which is quick, easy, and effective, with no chemicals involved. Contents stay sterile for 24 hours if the lid is unopened.

BPA Free



This sterilizer is made from BPA Free materials

Specifications

Weight and dimensions

- Dimensions: 290 x 160 x 150 (w x d x h) mm
- Weight: 1.5 kg

Country of origin

- Turkey

What is included

- Tongs: 1 pcs
- Newborn Pacifier: 1 pcs
- Classic Bottle 4oz: 1 pcs

Power

- Voltage: 220 V

Development stages

- Stages: 0 - 6 months

Compatibility

- Philips-AVENT range compatible

Material

- Polypropylene



Issue date 2013-04-04

Version: 5.1.1

© 2013 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

www.philips.com