



Gentle and accurate

Combining a digital and a soother thermometer.

Helping you to get a measurement of your baby's temperature, without upsetting your baby. The soother thermometer uses the same shield and teat as the AVENT Freeflow Soother for easy transition from soother to temperature measurement.

Comfort for your child

- Choice to adapt temperature measurement method to suit baby
- Soft, flexible tip digital thermometer
- Specially designed soother thermometer

Fast and reliable digital temperature readings

- Clear measurement displays and battery status indicators
- Battery low indicator

Convenient for you and your baby

- Includes AVENT Freeflow Soother

Specifications

Country of Origin

- Made in China: Thermometers, Yes
- Made in England: Soothers, Yes

Power

- Battery type: LR41
- Number of batteries: 1 per thermometer
- Operating time on battery: >100 hour(s)

Accessories

- User manual

Convenience

- Suitable for newborns
- Flexible tip
- Clear measurement display
- Automatic shut-off
- Battery status indicator

Technical specifications

- Operating temperature range: 32-42.9 °C
- Accuracy: 35-42 +/- 0.1 °C
- Measurement time: 10 secs ±3 secs in water (bath), 25 secs ±10 secs oral & rectal, 60 secs ±10 secs underarm



Issue date 2009-04-27

Version: 7.0.1

12 NC: 0000 000 00000

© 2009 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com