

# PHILIPS

## AVENT

Breast milk storage  
bags

Storage

Pre-sterilized

25 bags



SCF603/25



## Ultimate protection for your breast milk 6oz/180ml breast milk storage bags

Philips Avent breast milk storage bags provide safe and secure storage for your precious breast milk. They can be stored in the fridge or freezer and are pre-sterilized for immediate use.

### Pre-sterile and secure

- Freezer-safe with reinforced seams and dual layer bag
- Pre-sterilized bag with tamper proof seal
- Secure, leak-proof double zipper seal

### Easy to use

- Wide opening for easy filling and pouring
- The sturdy design of this bag enables it to be self standing
- Bags lay flat for easy storage

### Other benefits


- This bag is made out BPA free material



[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## Ultimate protection



Secure, leak-proof double zipper closure for safe storage of your breast milk

## Pre-sterilized, tamper proof




Tamper evident seal indicates the pre-sterilized bag is untouched before first use, for complete hygiene

## Strong and durable



Reinforced side seams and dual layered material provide extra reassurance that your precious breast milk is stored safely in the freezer or fridge

## Express, fill, store



The wide and sturdy opening of the bag ensures easy filling and pouring of your breast milk

## Self-standing bag



The sturdy design of this bag enables it to be self standing, and the wide opening ensures easy filling and pouring of breast milk

## Bags lay flat for easy storage

Bags lay flat for easy storage in freezer or fridge

## BPA free

This bag is made out BPA free material

# Specifications

## What is included

Bag: 25 pcs

## Capacity

Bag: 6oz/180ml

## Country of origin

China

## Material

BPA free\*

PET/PE

## Design

**Easy of use:** Wide, sturdy opening

**Secure:** Tamper evident seal

**Sturdy:** Self standing bag

**Ultimate protection:** Reinforced seams

## Features

**Leak proof:** Secure, double zipper

**Material:** Durable, dual layer bag

**Pre-sterilized**

## Dimensions

**Height:** 25,6 cm

**Width:** 9,9 cm



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.