

PHILIPS

AVENT

Digital baby
thermometer set

SCH540/02



Gentle and accurate

Combining a digital and a soother thermometer.

Helping you to get a measurement of your baby's temperature, without upsetting your baby. The soother thermometer uses the same shield and teat as the Avent Freeflow soother for easy transition from soother to temperature measurement.

Comfort for your child

- Choice to adapt temperature measurement method to suit baby
- Soft, flexible tip digital thermometer
- Specially designed pacifier thermometer

Fast and reliable digital temperature readings

- Clear measurement displays and battery status indicators
- Battery low indicator

Convenient for you and your baby

- Includes Avent Freeflow Pacifier

Specifications

Power

Operating time on battery: >100 hour(s)
Number of batteries: 1 per thermometer
Battery type: LR41

Accessories

User manual

Technical specifications

Operating temperature range: 32-42.9 °C
Measurement time: 10secs ±3secs in water (bath), 25secs ±10secs oral & rectal, 60secs ±10secs underarm
Accuracy: 35-42 +/- 0.1 °C

Convenience

Suitable for newborns
Flexible tip
Clear measurement display
Automatic shut-off
Battery status indicator

Country of origin

Made in China: Thermometers, Yes
Made in England: Soothers, Yes

Development stages

Stages: 0 - 6 months, 6 - 12 months

