

PHILIPS

AVENT

Digital baby
thermometer set

SCH540/00



Gentle and accurate

Combining a digital and a soother baby thermometer

Helping you to get a measurement of your baby's temperature, without upsetting it. The soother thermometer SCH540/00 uses the same shield and teat as the Avent Freeflow soother for easy transition from soother to temperature measurement.

Comfort for your child

- Choice to adapt temperature measurement method to suit baby
- Soft, flexible tip digital thermometer
- Specially designed pacifier thermometer

Fast and reliable digital temperature readings

- Clear measurement displays and battery status indicators
- Battery low indicator

Convenient for you and your baby

- Includes AVENT Freeflow Pacifier

Specifications

Country of origin

Made in China: Thermometers, Yes

Made in England: Soothers, Yes

Power

Battery type: LR41

Number of batteries: 1 per thermometer

Operating time on battery: >100 hour(s)

Accessories

User manual

Convenience

Suitable for newborns

Flexible tip

Clear measurement display

Automatic shut-off

Battery status indicator

Technical specifications

Operating temperature range: 32-42.9 °C

Accuracy: 35-42 +/- 0.1 °C

Measurement time: 10secs ±3secs in water (bath), 25secs ±10secs oral & rectal, 60secs ±10secs underarm

Development stages

Stages: 0 - 6 months, 6 - 12 months

