

# PHILIPS

## AVENT

### Microwave steam steriliser


Kills 99.9% of harmful germs

Sterilises in 2 minutes

Fits 4 Philips Avent bottles

Fits most microwaves

SCF281/03



## Fast microwave steam steriliser

For use at home or away

The Philips Avent microwave steam steriliser is a great choice for fast, effective sterilising at home or away. Up to four Philips Avent bottles or products can be sterilised at once, with 99.9% of germs killed in just 2 minutes.

#### **Fast and effective microwave steam sterilising**

- Microwave steam steriliser ready in just 2 minutes
- Kills 99.9% of germs and bacteria
- Contents stay sterile for up to 24 hours

#### **Fits various bottle types and other products**

- Sterilises standard neck and wide neck baby bottles
- Sterilises pumps, soothers, cutlery as well as baby bottles

#### **Compact, lightweight design**

- Lightweight design for sterile baby bottles on the go
- Fits most microwaves on the market
- Stay cool safety clips to keep lid secure



[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## Sterilises in just 2 minutes



The microwave steam steriliser can be used to sterilise baby bottles and other products in just 2 minutes, killing 99.9% of germs and bacteria. The exact length of the cycle depends on the wattage of your microwave. 2 minutes at 1200–1850 W, 4 minutes at 850–1100 W, 6 minutes at 500–800 W.

## Kills 99.9% of germs



The microwave steam steriliser is proven to kill 99.9% of germs and bacteria.

## Effective sterilisation



Contents remain sterile for up to 24 hours if the lid is unopened.

## Works across bottle types



The steriliser sterilises standard neck and wide neck baby bottles. Up to 4 Philips Avent baby bottles can be sterilised at once.

## Use for pumps, soothers too

The steriliser can be used to sterilise breast pumps, soothers, cutlery and other products as well as baby bottles.

## Lightweight design for travel



The steriliser's lightweight and compact design makes it easy to take along with you wherever you go — holidays, visits to relatives, etc. Now you can always have sterilised baby bottles and other products within easy reach.

## Fits most microwaves



The microwave steam steriliser has been designed to fit practically all microwaves on the market, keeping baby bottles and other products sterile wherever you are.

## Stay cool safety clips



The microwave steam steriliser has stay cool safety clips to keep the lid secure. This ensures you can remove the steriliser safely from the microwave, keeping baby bottles and other products sterile and preventing the possibility of burning.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Technical specifications

**Sterilisation time:** 2 min at 1200–1850 W, 4 min at 850–1100 W, 6 min at 500–800 W

**Water capacity:** 200 ml

## Weight and dimensions

**Weight:** 740 g

**Dimensions:** 166 (H), 280 (W), 280 (L) mm

## Country of origin

**Made in:** England

## What is included

**Microwave steam steriliser:** 1 pcs

**Tongs:** 1 pcs

## Development stages

**Stages:** 0–6 months

