

# PHILIPS

## AVENT

Digital baby  
thermometer set

SCH540/00



## Gentle and accurate

### Combining a digital and a soother baby thermometer

Helping you to get a measurement of your baby's temperature, without upsetting your baby. The soother thermometer SCH540/00 uses the same shield and teat as the Avent Freeflow soother for easy transition from soother to temperature measurement.

#### Comfort for your child

- Choice to adapt temperature measurement method to suit baby
- Soft, flexible tip digital thermometer
- Specially designed soother thermometer

#### Fast and reliable digital temperature readings

- Clear measurement displays and battery status indicators
- Battery low indicator

#### Convenient for you and your baby

- Includes Avent Freeflow Soother

# Specifications

## Country of origin

Made in China: Thermometers, Yes

Made in England: Soothers, Yes

## Power

Battery type: LR41

Number of batteries: 1 per thermometer

Operating time on battery: >100 hour(s)

## Accessories

User manual

## Convenience

Suitable for newborns

Flexible tip

Clear measurement display

Automatic shut-off

Battery status indicator

## Technical specifications

Operating temperature range: 32–42.9 °C

Accuracy: 35–42 +/- 0.1 °C

Measurement time: 10 secs ±3 secs in water (bath), 25 secs ±10 secs oral and rectal, 60 secs ±10 secs underarm

## Development stages

Stages: 0–6 months, 6–12 months

