

# PHILIPS

## AVENT

Digital baby  
thermometer set

SCH540/00



## Gentle and accurate

### Combining a digital and a soother baby thermometer

Helping you to get a measurement of your baby's temperature, without upsetting it. The soother thermometer SCH540/00 uses the same shield and teat as the Avent Freeflow soother for easy transition from soother to temperature measurement.

#### **Comfort for your child**

- Choice to adapt temperature measurement method to suit baby
- Soft, flexible tip digital thermometer
- Specially designed pacifier thermometer

#### **Fast and reliable digital temperature readings**

- Clear measurement displays and battery status indicators
- Battery low indicator

#### **Convenient for you and your baby**

- Includes Avent Freeflow Pacifier

# Specifications

## Country of origin

Made in China: Thermometers, Yes

Made in England: Soothers, Yes

## Power

Battery type: LR41

Number of batteries: 1 per thermometer

Operating time on battery: >100 hour(s)

## Accessories

User manual

## Convenience

Suitable for newborns

Flexible tip

Clear measurement display

Automatic shut-off

Battery status indicator

## Technical specifications

Operating temperature range: 32-42.9 °C

Accuracy: 35-42 +/- 0.1 °C

Measurement time: 10secs ±3secs in water (bath), 25secs ±10secs oral & rectal, 60secs ±10secs underarm

## Development stages

Stages: 0 - 6 months, 6 - 12 months

