

PHILIPS

AVENT

Baby Care set

All baby care essentials

Complete set

Boys and girls



SCH400/00



My first baby care set

All baby care essentials in one complete set

With our Baby care set SCH400 you can give the best care to your little one. The compact set includes an accurate, comfortable and quick-measuring digital thermometer, soft-tip nasal aspirator, finger toothbrush and nail and hair care set.

10 baby care essentials

- Digital thermometer with professional accuracy*
- Nasal Aspirator with soft and flexible tip
- Finger toothbrush
- Nail and hair care set

Compact pouch perfect for travelling

- Compact and structured set
- A perfect gift



asimpleswitch.com

Highlights

Digital thermometer



Quickly, comfortably and accurately measure your baby's temperature. The digital thermometer has professional accuracy* and the flexible tip provides extra comfort for you and your baby.

Compact and structured

The structured set provides room for your additional baby care products, whilst being compact and convenient to take with you. Ideal for travel and nursery.

Nasal aspirator



The nasal aspirator helps to open your baby's airways to help your baby feel and sleep better.

Finger toothbrush

Softly massage your baby's gums or teeth.

Nail and hair care set



Complete nail care set with rounded-tip scissors, nail clipper and three emery boards. And a hair care set with rounded-tip comb and soft hair brush.

Gift

A perfect gift for all (first-time) parents.

Specifications

Technical specifications

Battery: Battery 1.55 VDC (LR41)

Battery life: > 100 operating hours

Measurement range: 35°C to 42.9°C range

Accuracy: ±0.1°C for 35°C to 42°C range

Compliance: Conforms to MDD93/42/EEC

Logistic data

F-box dimensions (W x H x D): 117 x 174.5 x 51 mm

Development stages

Stages: All stages

Gender

For: Boy and girl



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

