

PHILIPS

AVENT

Manual breast pump with bottle

Easy for on-the-go use

Pump with soft massage cushion

Includes bottle and nipple

Travel cover and sealing disc



SCF330/20



More comfort, more milk, easy to take on-the-go*

Manual breast pump with massage cushion

Sit comfortably without leaning forward as our soft massage cushion gently stimulates your milk flow. The Avent manual breast pump has few parts, is easy to assemble, use & clean. Lightweight and compact, it's easy to express on-the-go.

Clinically proven less leaning forward*

- Relaxing position, thanks to its special design

Soft massage cushion with massaging petals

- Soft massage cushion gently stimulates milk flow

Easy to set up, express and clean

- Few separate parts and intuitive design
- Compact lightweight design

Easy to combine breast and bottle feeding

- Includes our Natural bottle and nipple for natural latch on



asimpleswitch.com

Highlights

Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, helps your milk to flow more easily.

Soft massage cushion



Our massage cushion has a soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let down.

Easy on-the-go



The breast pump is small and lightweight, which means it is easy to store and transport, making pumping on-the-go more discreet.

Easy set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

Includes our Natural bottle



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Specifications

Country of origin

England

Material

Breast pump: BPA free* (food contact parts only)

Bottle: BPA free*, Polypropylene

Nipple: BPA free*, Silicone

What is included

Breast pump body with handle: 1 pcs

Standard massage cushion: 1 pcs(19.5mm, a larger size cushion is available separately)

Natural bottle 4oz: 1 pcs

Extra soft newborn flow nipple: 1 pcs

Travel cover: 1 pcs

Sealing disc for milk storage: 1 pcs

Breast pad sample packs: 1 pack (2pcs)

Design

Breast pump design: Compact design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Development stages

Stages: 0 - 6 months

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation



[asimpleswitch.com](https://www.asimpleswitch.com)

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.



* Clinically proven comfort: Tests carried out on 110 mothers showed a significant preference for Philips Avent vs a leading competitor

* More Comfort: 73% out of 73 UK breastfeeding mothers agreed that this breast pump was more comfortable than their current pump (main market brands).

* Independent research has shown that there may be a link between stress levels and milk production. See www.philips.com/avent

* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011