

PHILIPS

AVENT

Silicone diaphragm

Restores suction to manual pump

Replacement diaphragm for Manual Pumps

SCF110/00



To bring vacuum in the breast pump.

Find matching products on the specifications tab

This silicone diaphragm creates (together with the stem) vacuum in your manual breast pump.

Ease of using Philips spare parts

- Easy renew your product with original Philips parts

Highlights

Easy renew of product parts

From time to time your product needs a facelift but it was never so easy like with Philips consumer replacement parts to renew your product. This all with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types :: SCF330

