

# PHILIPS

Infrared lamp



HP3614



## Relieves muscular pains and common colds

### Deep, focused warmth

Philips infrared lamp gives pain relief for muscles and joints. Its 100 Watt infrared warmth is comfortable and penetrates deep into skin, improving blood circulation thereby giving relief for focused areas of 20x30 cm.

#### **Relieves muscle pain**

- Deep penetrating Infrared warmth

#### **Focused treatment**

- 100 Watt infrared lamp

#### **Easy to use**

- Adjustable angle

#### **Easy to store**

- Cord storage

#### **Philips - the light expert**

- More than 15 million infrared lamps sold

# Highlights

## Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

## 100 Watt infrared lamp



The 100 Watt Philips Infrared lamp has been optimised for relief of pains in areas covering 20 x 30 cm, such as a shoulder, an elbow, a calf, or the neck.

## Adjustable angle



The infrared lamp can easily be positioned upto 40 degrees backwards.

## Cord storage



The infrared lamp has a cord storage holder at the back of the appliance. Simply wrap the cord around it and it can easily be stored again.

## Philips - the light expert

Philips is one of the global leading companies in light, with over 100 years of light experience. Globally Philips already sold more than 15 million of Infrared lamps.

# Specifications

## Easy storage

Cord storage

## Technical specifications

Power: 100 W

Voltage: 220/230 V

Frequency: 50 Hz

Cord length: 1.6 cm

Insulation: Class II (double isolation)

Type of lamps: R95, 100W

Lifetime of lamp: 750 sessions of 10 minutes

## Weight and dimensions

Product dimensions: 20.0x12.0x16.0 (HxWxD) cm

Product weight: 0.5 kg

F-box dimensions: 21.1x12.6x19.2 (HxWxD) cm

F-box weight: 1 kg

No. of F-boxes in A-box: 6 pcs

A-box dimensions: 22.4x26.2x58.8 (HxWxD) cm

A-box weight: 3.1 kg

Qty. on Euro pallet: 288 pcs

## Logistic data

CTV code: 884361401000

Country of origin: Germany

## Safety

IEC certified: In line with IEC 60335

## Easy to use

Adjustability: 0-40 degrees (backwards)

