

# PHILIPS

InfraPhil infrared  
lamp

150W



HP3616



## Relieves muscular pains

Deep, focused infrared warmth

Deep, focused warmth from infrared lamp

### **Relieves muscle pain**

- Deep penetrating Infrared warmth

### **Extra focused treatment**

- 150 Watt infrared lamp with extra focus

### **Easy to use**

- Adjustable angle
- On/off switch

### **Easy to store**

- Cord storage

# Highlights

## Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

## 150 Watt infrared lamp with extra focus



# 150W


The Philips Infrared lamp has been optimized for pain relief in areas covering 20 x 30 cm for example shoulder, elbow, calf, or neck. The concentric rings help to focus the beams on a specific area.

## Adjustable angle



The infrared lamp can easily be positioned upto 40 degrees backwards.

## On/off switch



Easy on/off switch at the back of the infrared lamp, so the appliance does not have to be unplugged after each treatment.

## Cord storage



The infrared lamp has a cord storage holder at the back of the appliance. Simply wrap the cord around it and it can easily be stored again.

# Specifications

## Easy to use

**Adjustability:** 0–40 degrees (backwards)  
**On/off switch**

## Easy storage

**Cord storage**

## Technical specifications

**Power:** 150 W

**Voltage:** 220/230 or 240/250V V

**Frequency:** 50 Hz

**Cord length:** 180 cm

**Insulation:** Class II (double isolation)

**Type of lamps:** PAR 38E, 150 W + prismatic rings for more focus

**Lifetime of lamp:** 750 sessions of 10 minutes

## Weight and dimensions

**Product dimensions:** 22.0x13.0x18.5 (HxWxD) cm

**Product weight:** 1 kg

**F-box dimensions:** 22.5x14.0x20.4 (HxWxD) cm

**F-box weight:** 1 kg

**No. of F-boxes in A-box:** 6 pcs

**A-box dimensions:** 24.0x28.8x62.4 (HxWxD) cm

**A-box weight:** 6.4 kg

**Qnt. on Euro pallet:** 192 pcs

## Logistic data

**CTV code:** 884361601000

**Country of origin:** Germany

## Safety

**IEC certified:** Complies to IEC 60335

## Medical appliance

**Medical Device Directive:** 2007/47/EC, MDD 93/42/EEC

