

PHILIPS

AVENT

Double Electric Breastfeeding set

Express, store and feed.

Includes 2 cushion sizes

Travel bag, handle, 4 bottles

25 Milk bags and 10 cups

SCF334/26



More milk in less time* Breastfeeding made easy

This set contains all you need to pump, store and feed your breast milk. Reduce expression time with double pumping, and sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

More milk in less time*

- Reduce expression time with double pumping

Effective expression with 3 customizable settings

- Choose the most effective and comfortable setting for you

Choose from a massaging and a power cushion

- Soft massage cushion with massaging petals
- Large power cushion for extra suction
- More comfortable pumping position due to unique design

Easy to set up, express and clean

- Few separate parts and intuitive design
- Includes our Natural bottle and nipple for natural latch on

Quiet pump for discreet expression, anytime

- Designed to let you express quietly



asimpleswitch.com

Highlights

More milk in less time



Reduce expression time with double pumping compared to single pumping.

Simple settings to choose from



When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortable for you.

Soft massage cushion



Our massage cushion has a soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let down.

Large power cushion



Choose the power cushion for the moments you need a stronger suction.

Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle or container, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, helps your milk to flow more easily.

Easy set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime. Closed system, no milk spillage in the tubes. Designed for hygienic expressing.

Includes our Natural bottle



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Quiet pump



Designed for discreet expression.



[asimpleswitch.com](https://www.asimpleswitch.com)

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Country of origin

England

Material

Breast pump: BPA free* (food contact parts only)

Bottle: BPA free*, Polypropylene

Nipple: BPA free*, Silicone

What is included

Breast pump body: 2 pcs

Base unit incl. tubing: 1 pcs

Manual handle: 1 pcs

Large massage cushion (25mm): 2 pcs (A smaller size cushion is available separately)

Power cushion (27mm): 2 pcs

Natural bottle 4oz: 4 pcs

Extra soft newborn flow nipple: 2 pcs

Travel cover: 2 pcs

Sealing disc for milk storage: 2 pcs

Milk storage cups: 10 pcs

Storage cups adaptors: 2 pcs

Milk storage bags: 25 pcs

Travel bag: 1 pcs

Breast pad sample packs: 1 pack (2pcs)

Design

Breast pump design: Compact design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Development stages

Stages: 0 – 6 months

Functions

More milk in less time: The perfect time saver

No leaning forward: Sit in a comfortable position

Settings: 3 Expression settings, 1 Stimulation mode

Soft massage cushion: Gentle stimulation



* More milk in less time: Reduce expression time with double pumping compared to single pumping.

* 60% of US breastfeeding mothers preferred the Philips Avent Comfort Double electric breast pump for comfort when expressing vs 29% of mothers who preferred the Medela in Style-Independent home placement test among 76 mothers in USA, August 2012. Medela is a registered trademark of Medela Holding AG

* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011