

# PHILIPS

## AVENT

Comfort Single  
electric breast pump

Natural

Includes 4 oz bottle



SCF332/01

## More comfort, more milk

### Breast pump with massage cushion

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

#### **Sit more comfortably with no need to lean forward**

- More comfortable pumping position due to unique design

#### **Choose your most comfortable setting**

- Features a gentle stimulation mode and 3 pumping settings

#### **Gently stimulates let-down and milk flow**

- Soft massage cushion with warm feel

#### **Easy to combine breast and bottle feeding**

- Includes our Natural bottle and teat for natural latch-on

#### **Other benefits**

- Compact lightweight design
- Intuitive assembly. Easy visual matching of parts
- Compatible with other Philips Avent feeding products
- Easy cleaning due to the small number of separate parts
- Easy operation at the touch of a button



[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means that you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

## Simple settings to choose from



When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow in the way that is most comfortable for you.

## Soft massage cushion



Our massage cushion has a new soft, velvety texture that gives the skin a warm feeling for comfortable, gentle stimulation of your milk

flow. The cushion is designed to gently mimic your baby's suckling to help stimulate let-down.

## Includes our Natural teat



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

## Compact design



The breast pump has a compact design, which makes it easy to hold and position on your breast. The small, lightweight base unit can easily be placed within comfortable reach for full control when pumping. For extra convenience in transport and storage, the tube simply wraps around the base unit. Easy to use on the go with batteries.

## Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. Your milk will never come into contact with the tubing and base unit. All parts are dishwasher proof, except the electrical parts.

## Easy operation

A single electric breast pump is ideal for mums who pump milk frequently and value being able to pump effortlessly, at the touch of a button.

## Intuitive assembly



Easy visual matching of parts for intuitive assembly.

## Fully compatible range



This breast pump can be used in combination with other feeding products in the Philips Avent range, including our Classic bottles and milk storage containers. Philips Avent also provides a range of breast care accessories to help you breastfeed for longer and enhance your comfort.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Country of origin

England

## Material

Breast pump: Polypropylene, BPA-free\*

Bottle: Polypropylene, BPA-free\*

Teat: Silicone, BPA-free\*

## What is included

Breast pump body: 1 pcs

Base unit incl. tubing: 1 pcs

Standard-size cushion: 1 pcs

Natural bottle 125 ml/4 oz: 1 pcs

Extra-soft newborn flow teat: 1 pcs

Sealing disc for milk storage: 1 pcs

Breast pad sample packs: 2 (2 Day pads and 2 Night pads) pcs

Spare diaphragm: 1 pcs

## Breast pump

Material: BPA-free\*

## Design

Breast pump design: Compact design

Bottle design: Ergonomic shape, Wide neck

## Ease of use

Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

## Development stages

Stages: 0–6 months

## Functions

No leaning forward: Sit in a comfortable position

Settings: 3 Expression settings, 1 Stimulation mode

Soft massage cushion: Gentle stimulation

