

PHILIPS

AVENT

Breast milk storage
bags

Storage

Pre-sterilized

25 bags



SCF603/25



Ultimate protection for your breast milk 6oz/180ml breast milk storage bags

Philips Avent breast milk storage bags provide safe and secure storage for your precious breast milk. They can be stored in the fridge or freezer and are pre-sterilized for immediate use.

Pre-sterile and secure

- Freezer-safe with reinforced seams and dual layer bag
- Pre-sterilized bag with tamper proof seal
- Secure, leak-proof double zipper seal

Easy to use

- Wide opening for easy filling and pouring
- The sturdy design of this bag enables it to be self standing
- Bags lay flat for easy storage

Other benefits


- This bag is made out BPA free material



asimpleswitch.com

Highlights

Ultimate protection




Secure, leak-proof double zipper closure for safe storage of your breast milk

Pre-sterilized, tamper proof



Tamper evident seal indicates the pre-sterilized bag is untouched before first use, for complete hygiene

Strong and durable



Reinforced side seams and dual layered material provide extra reassurance that your precious breast milk is stored safely in the freezer or fridge

Express, fill, store



The wide and sturdy opening of the bag ensures easy filling and pouring of your breast milk

Self-standing bag



The sturdy design of this bag enables it to be self standing, and the wide opening ensures easy filling and pouring of breast milk

Bags lay flat for easy storage

Bags lay flat for easy storage in freezer or fridge

BPA free

This bag is made out BPA free material

Specifications

What is included

Bag: 25 pcs

Capacity

Bag: 6oz/180ml

Country of origin

China

Material

BPA free*

PET/PE

Design

Easy of use: Wide, sturdy opening

Secure: Tamper evident seal

Sturdy: Self standing bag

Ultimate protection: Reinforced seams

Features

Leak proof: Secure, double zipper

Material: Durable, dual layer bag

Pre-sterilized

Dimensions

Height: 25,6 cm

Width: 9,9 cm



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.