

PHILIPS

AVENT

Handle for breast
pump

Natural

Comfort breast pump

Manual

CP9282



To make the breast pump diaphragm move.

Find matching products on the specifications tab

Pushing this handle up and down creates suction. You can control the suction by varying the pressure on the handle.

Ease of using Philips spare parts

- Easy renew your product with original Philips parts

Highlights

Easy renew of product parts

From time to time your product needs a facelift but it was never so easy like with Philips consumer replacement parts to renew your product. This all with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types :: SCF330

Fits product type :: SCF330/20, SCF330/40, SCF330/70, SCF330/71

