

















## For your child's safety and health

WARNING! • Before first and after every use, take all items apart, wash and rinse thoroughly. Sterilise using an AVENT Steam Steriliser or boil for five minutes. Wash your hands thoroughly and ensure surfaces are clean before contact with sterilised components. • DO NOT use abrasive cleaning agents or anti-bacterial cleaners. Excessive concentration of detergents may eventually cause plastic components to crack. Should this occur, replace immediately. • Dishwasher safe – food colourings may discolour components. 

Before each use, inspect all items. Throw away at the first sign of damage or weakness. For hygiene reasons, we recommend replacing teats/spouts after 3 months. • Always assemble teats when wet. • Keep teats in a dry, covered container, DO NOT leave in sunlight when not in use. Never use teats as a soother. • DO NOT warm contents in a microwave oven as this may cause uneven heating and could scald your baby. Always check food temperature before feeding. • Always use this product under adult supervision, DO NOT allow child to play with small parts or walk/run while using bottles or cups. • Continuous and prolonged sucking of fluids will cause tooth decay. Drinks other than milk and water, such as fruit juices and flavoured sugary drinks are not recommended. If they are used, they should be well diluted and only used for limited periods, not constantly sipped. 

MAGIC & SPORTSTER CUPS and spouts are not suitable for hot, fizzy or pulpy drinks. DO NOT use cups with spouts to mix infant formula as this may clog the non-spill valve and cause components to leak. Always ensure the valve is properly assembled. Magic cups should only be used as an aid to help children progress to using ordinary cups. . Conforms to EN 14350.

Spare teats and spouts with valves are available separately.

## AVENT® is here to help

Freephone: **0800 289064** (UK ONLY), **1800 509 021** (IRL ONLY), **1300 364 474** (AUSTRALIA NATIONAL), **02 9457 0127** (AUSTRALIA NSW ONLY). e-mail: england@cannonavent.com