



The Pacifier distributed in more than 2,000 hospitals nationwide.

Nipple specially designed to ensure natural development of teeth and gums. Pacifiers are suitable for sterilizing by means of a Philips AVENT Steam Sterilizer or boiling in water for at least 5 minutes. After sterilization please allow nipple to cool for at least 5 minutes and remove any trapped water before use.

IMPORTANT: AVENT pacifiers are available in sizes ranging from 0 to 18 months. For your child's safety, always follow age guidelines accordingly.

For your child's safety WARNING! Do Not Tie Pacifier Around Child's Neck as it Presents a Strangulation Danger.

In the event the pacifier becomes lodged in the mouth, DO NOT PANIC; it cannot be swallowed and is designed to cope with such an event. Remove from the mouth with care, as gently as possible. Always use this product under adult supervision. Inspect carefully before each use, especially when the child has teeth. Pull the pacifier in all directions. Throw away at the first signs of damage or weakness. We recommend replacing the pacifier after 4 weeks for safety and hygiene reasons. Do not leave a pacifier in direct sunlight or near a source of heat, or leave in disinfectant (sterilizing solution) for longer than recommended, as this may weaken the nipple. Never dip nipple in sweet substances or medication, your child may get tooth decay. Before first use, place in boiling water for 5 minutes (may cause slight color change), allow to cool and squeeze out any trapped water from the pacifier. This is to ensure hygiene. Clean before each use. The pacifier can be cleaned by washing with warm water; DO NOT use abrasive cleaning agents or antibacterial cleaners on the pacifiers. Excessive combinations of detergents may eventually cause components to crack. Should this occur, replace your pacifier immediately. Wash your hands thoroughly and ensure surfaces are clean before contact with sterilized components.

WARNING: Please check the nipple before use by pulling. Throw away at the first signs of damage or weakness. For safety and hygiene replace the pacifier after 4 weeks of use.

- 1 Hospital grade silicone.
- 2 One piece construction.
- 3 Highly durable.

1.800.54.AVENT
www.philips.com/AVENT

IN US MANUFACTURED FOR
PHILIPS CONSUMER LIFESTYLE,
A DIVISION OF PHILIPS ELECTRONICS
NORTH AMERICA CORPORATION,
P.O.BOX 10313, STAMFORD, CT 06904
MADE IN U.S.A.

Trademarks owned by the Philips Group. Soothie is a Registered Trademark of Respironics, Inc. and its affiliates.
© 2010 Koninklijke Philips Electronics N.V.
All Rights Reserved.



Pacifier - BPA-Free



4213 354 20511 SCF190/02

Conforms to US CPSC Standards