

## Philips Infrared lamps

The Philips InfraCare HP 3621 has been optimised for effective relief of pains in areas covering 20 x 30 cm, for example shoulder, elbow, calf, or neck. The InfraCare features innovative infrared halogen lamp technology. The special optics, filter and the powerful halogen lamp of 200 Watt have been optimised to provide a clear focused treatment as the infrared warmth is more equally divided over the treatment area.



## Our eco passport

Philips wants to make the world healthier and more sustainable. To make this a reality, we have to keep thinking differently; pushing the boundaries of what is currently believed possible. We continuously innovate and deliver products and technologies that are positive for people and the environment. At Philips we call this Green Innovation.

We know this is important to you too. So we have developed an eco passport which explains the environmental performance of our products in the Philips Green Focal Areas, developed through our eco design process since 1994, enabling you to make conscious, informed decisions. An eco passport will explain the Green Focal Areas considered and applied in the specific product.



## Energy

- Annual energy consumption <4kWh/year\*</li>
- · 200W infrared lamp
- 500 hours lamp lifetime



- Cardboard: contains >75% recycled fiber
- No use of polyvinylchloride (PVC) in packaging
- No use of expanded polystyrene (EPS) in packaging



## Substances

- The product is free of polyvinylchloride (PVC)\*\*
- The product is free of brominated flame retardants (BFR)\*\*

\*Energy consumption measured according to IEC 62301 and based on 50 times use per year, and 20 minutes usage per time. \*\*excluding power cord This passport applies to HP3621 and HP3643/01 range

Philips Green Product

2 years world-wide warranty



For service, spare and replaceable parts visit www.philips.com



Infrared Lamp



