

Philips Sonicare ProResults  
Compact sonic toothbrush  
heads

3-pack

*ProResults*

HX6023

## Focussed Clean

The smaller brush head size helps you easily maneuver the brush to focus on hard-to-reach areas.

### **Thorough clean feeling in your mouth**

- Angled, radial-trimmed bristles expands the area covered
- Contoured bristles fit the natural shape of your teeth
- Brush head design maximizes sonic motion

### **Easy to maneuver to focus on hard-to-reach areas**

- Small brush head size helps you focus in tight places

### **Designed to optimize performance**

- Reminder bristles ensure you are getting an effective clean

### **A brush head that fits multiple handles**


- Snap-on brush head

**PHILIPS**  
**sonicare**

sense and simplicity

## Highlights

### Snap-on brush head



Easy brush head placement and handle cleaning

### Angled radial-trimmed bristles




The curved surface area created by angled, radial-trimmed bristles helps maintain contact with the surfaces of the teeth and gums as the brush head pivots, covering a larger surface area without increasing the size of the brush head.

### Contoured bristles



Bristles are trimmed in a contoured pattern to fit the natural shape of teeth. The peaks and valleys along the length of the brush head are ergonomically designed to follow teeth topography and fit in between teeth so you're less likely to miss places when you brush.

### Maximize cleaning performance



Thoughtful design and precision engineering ensure that when this brush head is combined with the sonic motion and flexible brushing modes, you get a thorough clean.

### Reminder bristles

Reminder bristles let's you know when to replace the brush head. After three months of normal use bristles exhibit fatigue, and brush heads are less effective. Replace your brush heads every 3 months.

### Small brush head

The contoured bristle design on a smaller surface cleans as thoroughly as other Sonicare brush heads while allowing you to reach problem areas your dentist may have pointed out. This design helps you reach and focus on tight places where plaque and stains can hide.

## Specifications

### Ease of use

- Suitable for these models: FlexCare+, FlexCare, HealthyWhite, HydroClean, EasyClean, DiamondClean

### Technical specifications

- Operating time: For optimal results, a new brush head is recommended every three months.

### Weight and dimensions

- Dimensions brush head packaging: 18.3 h x 10.0 w x 4.8 d cm

### Items Included

- Hygienic travel cap: 3 pcs

