



Advanced Cleaning

Innovative technology delivers better results. Now there's a guaranteed way to better oral health. It's the new Philips Sonicare FlexCare that adapts to your oral care needs.

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Dynamic cleaning action drives fluid between teeth
- Removes up to 2x more plaque than a manual toothbrush
- Dynamic cleaning action drives fluid between teeth

Proven to improve oral health

- Naturally whiter teeth with patented sonic technology
- Improves gum health in only two weeks
- Helps whiten teeth
- Sensitive Mode: Gentle teeth and gum cleaning

Offers you a customized cleaning experience

- Easy-start program to acclimate to the Sonicare experience

Philips Sonicare FlexCare
Rechargeable sonic
toothbrush

3 modes

2 brush heads 1 soft travel case

FlexCare



HX6902\02

PHILIPS
sonicare

sense **and** simplicity

Highlights

Patented sonic technology



Sonicare's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Naturally whiter teeth



Dynamic Cleaning Action and direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

Improves gum health

This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks. Giving a superior cleaning in between teeth and removes significantly more plaque to a manual toothbrush

Smartimer



2-minute timer helps ensure dental professional recommended brushing time

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

Easy-start program

Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

Helps whiten teeth

Helps remove and reduce stains on your teeth for a brighter smile.

Removes up to 2x more plaque

Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

Gentle teeth and gum cleaning

Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums

Patented Sonic technology

Sonicare's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Specifications

Cleaning performance

- Performance: Removes up to 83% of plaque in hard-to-reach areas
- Health benefits: Helps improve gum health
- Whitening benefits: Helps remove tea, coffee and other stains
- Modes: 3 Modes: Clean, Sensitive, Massage and 2 Routines: Max Care and GoCare
- Brush heads: 1 ProResults standard brush head for an All-Around Clean and 1 ProResults compact brush head for a Precision Clean
- Speed: Up to 31,000 brush strokes per minute
- Timer: Smartimer and quadpacer

Ease of use

- Battery indicator: Illuminated icon indicates battery life, Light in power button blinks to indicate when to recharge
- Brush head system: Easy snap-on brush heads for optimal hygiene
- Brushing time: 28 2-minute brushings OR 2 weeks, 42 2-minute brushings OR 3 weeks
- Handle: Slim ergonomic design and rubber grip for easy handling.

Items included

- Additional brush head color code rings: 3 pcs
- Brush head color rings: 3 color rings for easy sharing
- Charger: 1 multi-voltage charger
- Charger cover: 1 charger cover with brush head holder
- Deluxe Charger with cord wrap: 1
- Deluxe soft travel case: 1
- FlexCare handle: 1

Technical specifications

- Battery: Rechargeable
- Operating time (full to empty): Fully charged should provide 20 2-min. brushings

Design and finishing

- Color: Clinical Green, White with sky blue grip and chrome front panel
- Chrome front panel

Service

- Warranty: 2-year limited warranty

Power

- Voltage: Multi-voltage charger

Features

- Deluxe recharge gauge: Light signal helps track battery charge
- Easy-Start: Gently increases power over 14 uses
- Quadpacer: Encourages attentive brushing of each quadrant of the mouth
- Rechargeable battery
- Smartimer: 2-minute timer helps ensure recommended brushing

