



Advanced Cleaning

Innovative technology delivers better results. Now there's a guaranteed way to better oral health. It's the new Philips Sonicare FlexCare that adapts to your oral care needs.

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

Offers you a customized cleaning experience

- Easy-start program to acclimate to the Sonicare experience

Proven to improve oral health

- Naturally whiter teeth with patented sonic technology
- Clinically proven safe and gentle
- Improves gum health in just two weeks

Philips Sonicare FlexCare
Rechargeable sonic
toothbrush

3 modes

2 brush heads 1 soft travel case

FlexCare



HX6932/75



PHILIPS
sonicare

sense and simplicity

Rechargeable sonic toothbrush
3 modes 2 brush heads 1 soft travel case

Specifications

Technical specifications

- Battery: Rechargeable

Design and finishing

- Color: White with sky blue grip and chrome front panel

Items included

- Charger: 1 multi-voltage charger with UV brush head sanitizer and 1 multi-voltage travel charger
- Brush head color rings: 3 color rings for easy sharing
- Sanitizer: 1 UV brush head sanitizer

Ease of use

- Brush head system: Easy snap-on brush heads for optimal hygiene
- Brushing time: 42 2-minute brushings OR 3 weeks
- Battery indicator: Illuminated icon indicates battery life
- Handle: Slim ergonomic design and rubber grip for easy handling.

Cleaning performance

- Speed: Up to 31,000 brush strokes per minute
- Brush heads: 1 ProResults standard brush head for an All-Around Clean and 1 ProResults compact brush head for a Precision Clean
- Performance: Removes up to 83% of plaque in hard-to-reach areas
- Health benefits: Helps improve gum health
- Whitening benefits: Helps remove tea, coffee and other stains
- Timer: Smarttimer and quadpacer
- Modes: 3 Modes: Clean, Sensitive, Massage and 2 Routines: Max Care and GoCare

Service

- Warranty: 2-year limited warranty

Power

- Voltage: Multi-voltage charger

HX6932/75

Highlights

Patented sonic technology

Sonicare's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Angled brush head neck

Uniquely angled brush head neck makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Naturally whiter teeth

Dynamic Cleaning Action and direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

Safe and gentle

Your Sonicare is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

Improves gum health

Optimal cleaning between teeth and along the gumline provides improved gum health in just two weeks. Giving a superior cleaning in between teeth and removes significantly more plaque to a manual toothbrush

Smartimer

Two-minute timer automatically tracks brushing time and turns the toothbrush off after two full minutes of brushing is complete.

Quadpacer

30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

Easy-start program

Gently increases power over first 14 uses to ease into the Sonicare experience



Issue date 2012-05-04

Version: 4.0.1

12 NC: 0000 000 00000

© 2012 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

www.philips.com