



# Complete Gum Care

Toothbrush for gentle cleaning of teeth and gum

Combining patented sonic technology and 5 brushing modes, FlexCare+ electric toothbrush is your solution for complete gum care. Now you can improve your gum health after just 2 weeks of use.

## Proven to improve oral health

- Naturally whiter teeth with patented sonic technology
- Clinically proven safe and gentle
- Improves gum health in only two weeks

## Provides a superior clean

- Dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

## Designed around you

- UV brush head Sanitizer - Cleaner brush, better oral health


## Helps you meet recommended brushing guidelines

- Quadpacer 30-second interval timer encourages even brushing
- Smartimer encourages 2-minute brushing time

Philips Sonicare FlexCare+  
Rechargeable sonic  
toothbrush with UV brush  
head Sanitizer

Dental Professional Trial

**FlexCare+**



HX6942/14



**PHILIPS**  
**sonicare**

sense and simplicity

## Highlights

### Patented sonic technology



Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

### Angled brush head neck

The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

### Naturally whiter teeth



The Dynamic Cleaning Action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.


### Safe and gentle

Your Philips Sonicare Electric Toothbrush is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

### Improves gum health

This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks. Giving a superior cleaning in between teeth and removes significantly more plaque to a manual toothbrush

### Quadpacer



Quadpacer 30-second interval timer encourages even brushing

### Smartimer



2-minute timer helps ensure dental professional recommended brushing time

### UV brush head Sanitizer

All-in-one sanitizing, charging and storage. The UV brush head Sanitizer helps keep your brush head free of bacteria. It kills up to 99% of bacteria and viruses.\*

## Specifications

### Cleaning performance

- Performance: Removes up to 83% of plaque in hard-to-reach areas
- Health benefits: Improves gum health in only two weeks
- Whitening benefits: Helps remove tea, coffee and other stains
- Modes: 5 Modes: Clean, Sensitive, Massage, Refresh and Gum Care
- Speed: Up to 31,000 brush strokes per minute
- Timer: Smartimer and quadpacer

### Ease of use

- Battery indicator: Illuminated icon indicates battery life
- Brush head system: Easy snap-on brush heads for optimal hygiene
- Brushing time: Up to 3 weeks
- Display: Illuminated display indicates selected mode, becomes invisible when power is off.
- Handle: Slim ergonomic design and rubber grip for easy handling.

### Items included

- Brush head color rings: 3 color rings for easy sharing
- Charger: 1 multi-voltage charger with UV brush head sanitizer and 1 multi-voltage travel charger
- Sanitizer: 1 UV brush head sanitizer
- Brush heads: 1 ProResults compact, 1 ProResults standard

### Technical specifications

- Battery: Rechargeable

### Design and finishing

- Color: Soft Sky Blue, White with sky blue grip and illuminated display

### Service

- Warranty: 2-year limited warranty

### Power

- Voltage: Multi-voltage charger



Issue date 2013-04-04

Version: 2.0.1

12 NC: 0000 000 00000

© 2013 Koninklijke Philips Electronics N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

[www.philips.com](http://www.philips.com)

\* E coli, S. mutans and Herpes simplex