



Philips Sonicare Essence  
Rechargeable sonic  
toothbrush

**1 mode**

2 brush heads

**HX5352**

## Better plaque removal

Sonicare's unique dynamic cleaning action gently and effectively reaches deep between teeth and along the gumline.

### **Guides you to follow dental recommendations**

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

### **Provides a superior clean**

- Dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

### **Proven to improve oral health**

- Clinically proven safe and gentle

### **Offers a customized brushing experience**

- Easy-start program to acclimate to the Sonicare experience

**PHILIPS**  
**sonicare**

sense and simplicity

Rechargeable sonic toothbrush  
1 mode 2 brush heads

# Specifications

## Technical specifications

- Battery: Rechargeable

## Design and finishing

- Color: White

## Items included

- Charger: 1 charger
- Charger cover: 1 charger cover with brush head holder

## Ease of use

- Brush head system: Easy screw-on brush heads
- Brushing time: 28 2-minute brushings OR 2 weeks
- Battery indicator: Light in power button blinks to indicate when to recharge

## Cleaning performance

- Speed: Up to 31,000 brush strokes per minute
- Brush heads: 2 e Series standard brush heads
- Performance: Better plaque removal than a manual toothbrush
- Health benefits: Helps improve gum health
- Whitening benefits: Helps remove tea, coffee and other stains
- Timer: Smartimer and quadpacer

## Service

- Warranty: 2-year limited warranty

## Power

- Voltage: 110-120 V, 50-60 Hz

HX5352/55

# Highlights

## Patented sonic technology



Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.


## Angled brush head neck

The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Safe and gentle

Your Philips Sonicare Electric Toothbrush is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

## Smartimer



2-minute timer on this toothbrush helps ensure dental professional recommended brushing time

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

## Easy-start program

Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience



Issue date 2013-04-04

Version: 4.1.1

© 2013 Koninklijke Philips Electronics N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

[www.philips.com](http://www.philips.com)