

Philips Sonicare
Battery Sonic Toothbrush



HX3351



Longer brushing *1



Safe on braces and brackets *4

Proven to deliver better brushing compliance for teenagers

Encourages better brushing habits

Putting your child on the right track for oral care at an early age reinforces positive, long-lasting healthcare habits throughout life, even when you are not there.

Gives your child healthier teeth and gums

- Longer brushing *1
- Removes significantly more plaque *2
- Positive experience *3
- Safe on braces and brackets *4

PHILIPS
sonicare
the sonic toothbrush

PHILIPS
sense and simplicity

Specifications

Technical Specifications

- **Speed:** 31,000 brush strokes per minute
- **Battery type:** 1.5 V AA alkaline batteries (2 included)

Accessories

- **Brush head:** One standard e-series contoured brush head

Design and finishing

- **Material handle:** ABS
- **Material brush head:** PA, PCTA, TPE with nylon bristles
- **Colour:** Blue

Weight and dimensions

- **F-box dimensions:** 4.3 x 2.3 x 9.4 in
- **F-box weight:** 0.57 lb
- **A-box dimensions:** 9.8 x 6.1 x 4.8 in
- **A-box weight:** 1.96 lb
- **A-box number of pieces:** 3 pcs

References

- *1,2,3,4: Data on file. 2005

Product highlights

Longer brushing *1

Tweens and teens brush longer using Xtreme than a manual toothbrush.

Removes significantly more plaque *2

Clinically proven to remove significantly more plaque than a manual toothbrush.

Positive experience *3

97% of tweens and teens rated Xtreme favourably in overall liking.



Issue date 2009-02-16

Version: 2.2

© 2009 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com