

Philips Sonicare EasyClean  
Sonic electric toothbrush

2 Modes

*EasyClean*

HX6582

## Better plaque removal

Sonicare's unique dynamic cleaning action gently and effectively reaches deep between teeth and along the gumline.

### **Guides you to follow dental recommendations**

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

### **Provides a superior clean**

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

### **Proven to improve oral health**

- Clinically proven safe and gentle

### **Offers a customized brushing experience**

- Easy-start program to acclimate to the Sonicare experience



[asimpleswitch.com](http://asimpleswitch.com)

# PHILIPS

**sonicare**

## Highlights

### Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

### Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

### Safe and gentle




Your Philips Sonicare Electric Toothbrush is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

### Smartimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

### Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

### Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

## Specifications

### Technical specifications

- Operating time (full to empty): Up to 2 weeks

### Ease of use

- Brushing time: Up to 2 weeks



Issue date 2016-10-29

Version: 1.0.1

© 2016 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)