



Philips Sonicare FlexCare+  
Sonic electric toothbrush

**5 modes**

2 brush heads  
UV brush head sanitizer

**FlexCare+**



**HX6972/09**

## Complete gum care

Toothbrush for gentle cleaning of teeth and gum

Combining our sonic technology and 5 brushing modes, FlexCare+ electric toothbrush is your solution for improved oral health. Plus you can improve your gum health after just 2 weeks of use.

### **Offers a customized brushing experience**

- Easy-start program to acclimate to the Sonicare experience
- Five different modes for a better cleaning experience

### **Provides a superior clean**

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

### **Guides you to follow dental recommendations**

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

### **Proven to improve oral health**

- Naturally whiter teeth with patented sonic technology
- Clinically proven safe and gentle
- Improves gum health in only two weeks



# PHILIPS

**sonicare**

# Highlights

## Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

## Safe and gentle



Your Philips Sonicare Electric Toothbrush is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

## Improves gum health



This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks.

## Smartimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

## Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

## Five different cleaning modes



The five modes include a Clean mode for an outstanding clean in two minutes. A Sensitive mode for gentle yet thorough cleaning of the gums. A Refresh mode for those quick one minute touch-ups. A Massage mode to stimulate the gum and a Gum Care mode that helps improve gum health.

# Specifications

## Modes

- Clean: For exceptional everyday clean
- Gum Care: Gently massages gums
- Refresh: 1 min mode for quick touch up
- Sensitive: Gentle teeth and gum cleaning
- Massage: Invigorating massage

## Items included

- Handles: 1 FlexCare+
- Brush heads: 1 ProResults standard, 1 ProResults compact
- Charger: 1
- UV sanitizer
- Travel case: 1

## Design and finishing

- Color: White with soft sky blue

## Cleaning performance

- Speed: Up to 62000 brush movement/min
- Performance: Removes up to 83% of plaque in hard-to-reach areas
- Health benefits: Helps improve gum health

- Whitening benefits: Helps naturally remove stains
- Timer: Smartimer and Quadpacer
- Pressure sensor: Vibrates handle to alert user

## Ease of use

- Brush head system: Easy click-on brush heads
- Brushing time: Up to 3 weeks\*\*
- Battery indicator: Light shows battery status
- Handle: Slim ergonomic design, Rubber grip for easy handling
- Display: Illuminated display

## Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Operating time (full to empty): Up to 3 weeks\*\*

## Power

- Voltage: 110-220 V

## Service

- Warranty: 2-year limited warranty



Issue date 2018-03-11

Version: 6.0.1

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* E coli, S. mutans and Herpes simplex

\*\*based on two periods of two-minute brushings per day, on clean mode