

# PHILIPS

## sonicare

Sonic electric  
toothbrush - Dispense

HealthyWhite

2 modes

2 brush heads



HX6782/34



## Whiter, healthier teeth

Everybody loves a bright white smile. Let the natural whiteness of your patients' teeth shine through. HealthyWhite is proven to remove everyday stains in just two weeks, with the regular use of the Clean & White Mode.

### Proven to improve oral health

- Naturally whiter teeth with patented sonic technology

### Offers a customized brushing experience

- Clean and White Mode: proven to remove stains

### Designed around you

- EasyStart

### Guides you to follow dental recommendations

- Quadpacer
- Smartimer

### Provides a superior clean

- Safe on orthodontics, dental work and veneers
- Angled brush head neck for better reach to back teeth
- Engineered to maximize sonic motion



asimpleswitch.com

# Highlights

## Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

## Safe and gentle



The gentle motion of this Philips Sonicare power toothbrush is safe for dental work, including orthodontics, implants and veneers, so you can brush without worry.

## Clean and White Mode



2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.\*

## Sonic technology



Philips Sonicare brush heads are vital to our core technology of high-frequency, high-amplitude brush movements performing over 31,000 brush strokes per minute. This dynamic cleaning action drives fluid deep between teeth and along the gumline, for a thorough yet gentle cleaning.

## Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Quadpacer




Thirty-second interval timer encourages thorough brushing of each quadrant of the mouth.

## Smartimer



Two-minute timer helps ensure dental professional recommended brushing time.

## EasyStart



Gently increases power over the first 14 uses to ease patients into the Sonicare brushing experience.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Technical specifications

**Battery:** Rechargeable

## Modes

**Clean:** For exceptional everyday clean

**White:** Removes surface stains

## Design and finishing

**Color:** Glacier Green

## Items included

**Brush heads:** 2 ProResults standard

**Charger:** 1

**Handles:** 2 HealthyWhite

**Travel case:** 1

## Cleaning performance

**Health benefits:** Helps improve gum health

**Performance:** Better plaque removal than a manual toothbrush

**Speed:** Up to 62000 brush movement/min

**Timer:** SmartTimer and Quadpacer

**Whitening benefits:** Helps naturally remove stains

## Ease of use

**Battery indicator:** Light shows battery status

**Brushing time:** Up to 3 weeks\*\*

**Handle:** Slim ergonomic design

## Power

**Voltage:** Multi-voltage charger

## Service

**Warranty:** 2-year limited warranty

