





## Advanced plaque removal made simple

The switch from a manual toothbrush made easy

Philips Sonicare Xtreme is the battery-powered sonic toothbrush that is clinically proven to deliver superior plaque removal versus a manual toothbrush

## Provides a superior clean

- Soft bristles
- · Patented Sonicare toothbrush technology

## Proven to improve oral health

· Clinically proven safe and gentle

## Guides you to follow dental recommendations

Smartimer

# Highlights

### Safe and gentle



Your Philips Sonicare Electric Toothbrush is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

### **Smartimer**



Two-minute timer helps ensure dental professional recommended brushing time.

#### Soft bristles

Philips Sonicare e-Series brush heads have extra gentle bristles to ensure deep yet gentle cleaning.

### Sonic technology



Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

# Specifications

### Modes

Clean: For exceptional everyday clean

### **Items included**

Handles: 1 Essence

Brush heads: 1 e-Series standard

Charger: 1

### **Design and finishing**

Color: White

### **Cleaning performance**

Speed: Up to 62000 brush movement/min Performance: Better plaque removal than a

manual toothbrush

Whitening benefits: Helps naturally remove stains

Timer: SmarTimer

## **Ease of use**Brush head system: Easy screw-on brush

heads

Brushing time: Up to 2 weeks Handle: Ergonomic design

Battery indicator: Light shows battery status

### Power

Voltage: 110-220 V

### **Service**

Warranty: 2-year limited warranty

## **Technical specifications**

Battery: Two AA batteries included

Battery type: NiMH

Operating time (full to empty): Up to 2 weeks

© 2019 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2019-07-17 Version: 3.1.1

12 NC: 0000 000 00000 EAN: 00 07502 08033 51 www.philips.com

