

# PHILIPS

## sonicare

Sonic electric  
toothbrush

Elite+

1 mode

1 brush head



HX3216/13



### Clinically proven to provide superior cleaning \*

### Up to 3x more plaque removal\*

Designed to be the perfect toothbrush to transition from your manual toothbrush. Until now, it feels like you've never brushed your teeth before.

#### **Innovative technology**

- Philips Sonicare's advanced sonic technology

#### **Proven to improve oral health**

- Up to 3x more plaque removal than a manual toothbrush

#### **Designed around you**

- Timers encourage thorough 2 minutes of brushing
- A great way to get used to your new toothbrush

#### **Optimize your brushing**

- Battery life
- Click-on brush head system



[asimpleswitch.com](http://asimpleswitch.com)


# Highlights

## Philips Sonicare technology




Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

## Say goodbye to plaque



Densely-packed, high-quality bristles remove up to 3x more plaque than a manual toothbrush. And the contoured bristles fit the shape of your teeth, so you can easily clean those hard-to-reach areas.


## QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our

Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

## Easy-start program



Starting something new can take getting used to. That's why our easy-start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

## Battery life



Lasts up to 14 days between charges.

## Click-on brush head system



This brush head clicks on and off your brush handle for a secure fit and easy maintenance and cleaning. It fits Philips Sonicare toothbrush PowerUp Battery, Essence+/Elite+/CleanCare+ handles.

# Specifications

## Modes

**Clean:** For exceptional everyday clean

## Items included

Handle: 1 Elite+

Brush heads: 1 C1 ProResults standard

Charger: 1

## Design and finishing

Color: Mid-Blue

## Cleaning performance

Timer: Smartimer and Quadpacer

Performance: Removes up to 3x more plaque\*

## Ease of use

Handle: Ergonomic design

Brush head system: Easy click-on brush heads

## Technical specifications

Battery: Rechargeable

Battery type: NiMH

Operating time (full to empty): Up to 14 days

## Power

Voltage: 110-220 V

## Service

Warranty: 2-year limited warranty



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.