

# PHILIPS

## sonicare

Battery Sonicare  
toothbrush

PowerUp

1 mode

1 brush head

HX3631/02



## Simple, Gentle, Superior

More strokes in a day than your manual in a month

Designed to be as easy to use as your manual toothbrush, while delivering superior results and a better experience. With more brush strokes in a day than your manual in a month, PowerUp removes more plaque, helps reduce cavities and improve gum health

### Gentle

- Clinically proven safe and gentle

### Simple

- Familiar brushing motion

### Superior

- Sonic Technology
- Removes more plaque than a manual toothbrush
- 2 minute timer helps ensure recommended brushing time
- Helps improve gum health
- Helps reduce cavities


# Highlights

## Familiar brushing motion



Familiar brushing motion and brush head shape of a manual toothbrush, but with the cleaning power of over 15,000 strokes per minute

## Helps improve gum health



Effectively helps improve gum health via a gentle cleaning action

## Helps reduce cavities



Brushing twice a day every day helps reduce cavities

## PowerUp Sonic Technology



More strokes in 1 day than your manual in a month. Over 15,000 strokes per minute

## Removes more plaque



Removes more plaque than a manual toothbrush.

## Safe and Gentle



Your toothbrush is safe to use on all dental work

## Smart timer



Timer helps ensure dental professional recommended brushing time.

# Specifications

## Modes

**Clean:** For exceptional everyday clean

## Items included

**Handles:** 1 PowerUp

**Brush heads:** 1 PowerUp brush head

**Batteries:** 2 AA

## Design and finishing

**Color:** Blue

## Cleaning performance

**Speed:** Up to 15,000 brush strokes/min

**Performance:** Helps remove plaque

**Health benefits:** Helps improve gum health

**Whitening benefits:** Helps whiten teeth

**Timer:** Smartimer

## Ease of use

**Brush head system:** PowerUp brush heads

**Battery Life:** 120 2-minute brushings

**Handle:** Ergonomic design, Rubber grip for easy handling

## Technical specifications

**Battery:** 2 AA Battery

**Battery type:** Alkaline

## Service

**Warranty:** 1-year limited warranty

