

# PHILIPS

## sonicare

Rechargeable  
toothbrush

Sensiflex

1 mode



HX1620



## Naturally whiter teeth

### Toothbrush for whiter teeth

Your smile says it all. What better way to get naturally whiter teeth than with Philips Sensiflex 1620 electric toothbrush? The Dual Cleaning Action brush head delivers a clinically better clean compared to a manual toothbrush.

#### Double Cleaning Action

- Naturally whiter teeth
- Cleans visible teeth surfaces
- Targets hard-to-reach plaque in between teeth

#### Gentle on gums, tough on plaque

- Gum Protection system
- Adjust to maintain optimal brushing pressure

# Highlights

## Double Cleaning Action



Cleans visible teeth surfaces while the active tip targets hard-to-reach plaque in between teeth.

## Gum Protection system



Philips Sonicare's patented Gum Protection system adjusts to maintain optimal brushing pressure

# Specifications

## Cleaning performance

Brush heads: 1

## Items included

Charger base with brush head storage pins

## Technical specifications

**Battery:** Rechargeable

**Charging time:** 16 hour(s)

**Cord length:** 1.16 m

**Dual brush head:** Brush head combines to create 15,000 strokes/min.

**Operating time (full to empty):** 25 minute(s)

**Power consumption:** 2.9 W

