

PHILIPS

sonicare

Sonic electric
toothbrush - Dispense

Essence+

1 mode

1 brush head



HX3281/32



Removes up to 3x more plaque* for a superior clean

*than a manual toothbrush

With its one button mode and EasyStart feature, this Philips Sonicare electric toothbrush is the simple way to help patients find a gentler power toothbrush experience. It removes up to 3x more plaque than a manual toothbrush.

Provides a superior clean

- Optimized to deliver superior plaque removal*
- Patented Sonicare toothbrush technology

Designed around you

- EasyStart



asimpleswitch.com

Highlights

Sonic technology




Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

Superior plaque removal



This brush head features densely-packed, high-quality tufts to effectively remove plaque. Also features specially curved power tufts to easily reach teeth at the back of the mouth.

EasyStart



Gently increases power over the first 14 uses to ease patients into the Sonicare brushing experience.

Specifications

Modes

Clean: For exceptional everyday clean

Items included

Brush heads: 1 ProResults plaque control

Charger: 1

Handles: 1 Essence+

Design and finishing

Color: Fresh green

Cleaning performance

Health benefits: Helps improve gum health

Performance: Removes up to 3x more plaque*

Speed: Up to 31000 brush movement/min

Timer: SmarTimer

Whitening benefits: Helps maintain white teeth

Ease of use

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Brushing time: Up to 10 days

Handle: Ergonomic design, Rubber grip for easy handling

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

Technical specifications

Battery: Rechargeable

Battery type: NiMH

Operating time (full to empty): Up to 10 days



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.