

PHILIPS

sonicare

Sonic electric
toothbrush

FlexCare Platinum

3 modes, 3 intensities

2 brush heads

With pressure sensor



HX9112/52



Deep, gentle cleaning

Maintains gum health & 10x more plaque removal than a manual toothbrush

With FlexCare Platinum, even the deepest clean is gentle enough to protect your teeth and gums. Thanks to the AdaptiveClean brush head and intuitive pressure sensor, you get maximum plaque removal and a gentle experience however you brush.

Proven to improve oral health

- Remove up to 10x more plaque with AdaptiveClean (versus a manual toothbrush)
- Safe and gentle for your oral health needs

Optimize your brushing

- Choose from 3 modes and 3 intensity settings

Provides a superior clean

- AdaptiveClean brush head flexes along your gum line
- Deep Clean mode tackles your trouble spots
- Our unique technology brings you a powerful yet gentle clean
- Angled neck to help tackle cleaning in tough areas

Designed around you


- Alerts you if you're pressing too hard
- Designed to fit your lifestyle
- Easy start program builds your Philips Sonicare routine
- Timers encourage thorough 2 minutes of brushing



asimpleswitch.com

Highlights

Remove up to 10x more plaque versus a manual toothbrush



Click on the AdaptiveClean brush head to experience our deepest ever clean. You get 4x more surface contact** to remove up to 10x more plaque in hard-to-reach areas*, and a truly deep clean along the gum line and between teeth.*vs a manual toothbrush**as compared to DiamondClean brush head

Up to 2x healthier teeth




Click on our AdaptiveClean brush head for a deep, gentle clean that helps maintain gum health. Soft rubber sides flex along every contour of your teeth and gums to absorb any excess pressure. Bristles can then reach deep between teeth giving you a comfortable, effective deep clean.

Extra deep cleaning



Our Deep Clean mode means you can really focus on specific areas, paying them extra attention. We've combined our unique brush head motion with an extended 3 minutes of brushing time so your teeth get a truly deep clean.

3 modes, 3 intensity settings



With FlexCare Platinum, you experience a truly deep clean. Our 3 intensity settings let you boost your clean, while 3 modes meet your specific needs: Clean mode – for superior daily cleaning, White – to remove surface stains, and Deep Clean for invigorating freshness.

Too much pressure?



You might not notice if you're brushing too hard, but your FlexCare Platinum will. If you use too much pressure, your handle will gently pulsate. This is a reminder to ease off the pressure, and let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

Long battery life



Our premium travel case lets you store your toothbrush hygienically, while our compact travel charger keeps you topped up when you're on the go. You'll enjoy two weeks of regular use from a single full charge, so being on the move means you needn't miss out on that fresh feeling.

Ease into your new routine



Starting something new can take some getting used to. Our easy-start program gives you the option of a gradual, gentle increase in brushing power with the first 14 uses with your new toothbrush.

QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Modes

Clean: For exceptional everyday clean

White: Removes surface stains

Deep Clean: For an invigorating deep clean

3 intensities: High, Low, Medium

Items included

Handle: 1 FlexCare Platinum

Brush heads: 1 AdaptiveClean, 1 InterCare standard

Travel case: 1

Charger: 1

Design and finishing

Color: Slate gray

Cleaning performance

Speed: Up to 62000 brush movement/min

Performance: Removes up to 10x more plaque versus a manual toothbrush

Health benefits: Maintains gum health

Whitening benefits: Helps naturally remove stains

Timer: SmarTimer and Quadpacer

Pressure sensor: Vibrates handle to alert user

Ease of use

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Display: Illuminated display

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): Up to 3 weeks

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

