

# PHILIPS

## sonicare

Sonic electric  
toothbrush

HealthyWhite

3 Modes

2 brush heads



HX6732/45



## Whiter, healthier teeth

### The ultimate toothbrush for whiter teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. Sonicare HealthyWhite electric toothbrush is proven to remove everyday stains in just two weeks, with the regular use of the Clean and White Mode.

#### Proven to improve oral health

- Naturally whiter teeth with patented sonic technology
- Clinically proven safe and gentle

#### Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

#### Offers a customised brushing experience

- Easy-start programme to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains
- Clean Mode: Standard 2-minute mode
- Sensitive Mode: Gentle teeth and gum cleaning

#### Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth



asimpleswitch.com

# Highlights

## Sonic technology




The Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.

## Angled brush head neck



The uniquely angled neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has been proven to remove everyday stains, to help keep teeth naturally whiter.

## Safe and gentle



Your Philips Sonicare Electric Toothbrush is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

## Smartimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

## Easy-start programme



Gently increases power over first 14 uses to ease you into the Philips Sonicare electric toothbrush experience

## Clean and White Mode



2 minutes of Clean mode with an additional 30 seconds of White mode to focus on your visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.\*

## Clean Mode

For exceptional everyday cleaning. Helps maintain whitening treatments.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Modes

**Clean:** For exceptional everyday cleaning

**Sensitive:** Gentle teeth and gum cleaning

**Clean and White:** 2 min Clean + 30 sec White

## Items included

**Handles:** 1 HealthyWhite

**Brush heads:** 1 ProResults standard, 1  
DiamondClean standard

**Charger:** 1

**Travel case:** 1

## Design and finishing

**Colour:** Frost White

## Cleaning performance

**Speed:** Up to 31,000 brush strokes/min

**Performance:** Removes up to 90% of stains\*

**Health benefits:** Helps improve gum health

**Whitening benefits:** Whitens teeth by up to 2 shades

**Timer:** Smartimer and Quadpacer

## Ease of use

**Brush head system:** Easy click-on brush heads

**Brushing time:** Up to 3 weeks\*\*

**Battery indicator:** Light shows battery status

**Handle:** Slim ergonomic design, Rubber grip for easy handling

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Operating time (full to empty):** Up to 3 weeks\*\*

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

