

PHILIPS

sonicare

Sonic electric
toothbrush - Dispense

DiamondClean

5 modes

4 brush heads

Glass charger, travel case
with Deep Clean mode



HX9383/76



Ultimate cleaning, superior results

Superior plaque removal plus powerful whitening

Only Philips Sonicare DiamondClean provides the performance your patients need and the features they want. DiamondClean removes up to 10x more plaque than a manual brush and eliminates surface stains to whiten smiles in just one week.

Provides a superior clean

- DiamondClean toothbrush head for Sonicare's best whitening
- AdaptiveClean brush head
- New Deep Clean: Effortlessly tackle trouble spots
- Philips Sonicare toothbrush technology

Guides you to follow dental recommendations

- 5 modes: Clean, White, Sensitive, Gum Care, Deep Clean

Designed around you

- Illuminated display
- Three-week battery life
- Dual charging system


Highlights

DiamondClean brush head




DiamondClean toothbrush head, Philips Sonicare's best whitening brush head, has diamond shaped, medium stiffness bristles to effectively yet gently scrape away plaque. This Philips Sonicare electric toothbrush provides an exceptional cleaning and whiter teeth compared to a manual toothbrush.

AdaptiveClean brush head




Soft, flexible sides made of rubber allow the bristles to conform to the unique shape of patients' teeth and gums to provide up to 4x more surface contact*. The brush head's flexibility helps the bristles gently track along the gumline and absorb excessive brushing pressure. Enhanced sweeping motion delivers a unique in-mouth feel and superior cleaning along the gumline and in hard-to-reach areas*. Proven to remove 10x more plaque than a manual toothbrush in hard-to-reach areas*.

Five brushing modes



Includes the stain-removing White mode, clinically proven to whiten teeth in as little as one week.

Deep Clean mode




Our all-new Deep Clean mode delivers an exceptional clean in hard-to-reach places. The unique motion of the brush head in this mode plus an extended three-minute brushing time makes it the ideal choice when you want to pay extra attention to trouble spots.

Charging system




Innovative glass charger and station for home use, plus portable charging case with laptop-ready USB and wall plug for travel.

Illuminated display



When on, unique backlight indicates selected mode. When off, display becomes invisible for a clean look.

Rechargeable battery



When fully charged, DiamondClean has enough battery life for three weeks of regular use.

Sonic technology



Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

Specifications

Modes

Clean: For exceptional everyday clean

White: Removes surface stains

Sensitive: Gentle teeth and gum cleaning

Gum Care: Gently massages gums

Deep Clean: For an invigorating deep clean

Items included

Handles: 1 DiamondClean

Brush heads: 1 InterCare standard, 2

DiamondClean standard, 1 AdaptiveClean

Glass charger: 1

Travel case: USB travel charger

Design and finishing

Color: Amethyst

Cleaning performance

Speed: Up to 62000 brush movement/min

Health benefits: Improves gum health in only two weeks

Performance: Removes up to 10x more plaque*

Whitening benefits: Whitens teeth up to 2 shades

Timer: SmarTimer and Quadpacer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks**

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Display: Illuminated display

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): 3 week battery life

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

